NAUSICAA

MEDICAL

www.nausicaa-medical.com

FRANCAIS

ENGLISH

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FRANCAIS

ENGLISH

PATIENT LIFTS : NAUSIFLY 4 COMPACT AND STANDARD

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EN

USER MANUAL

HOIST: **NAUSIFLY 4 COMPACT AND STANDARD WITH ECP MECHANICAL OR ELECTRICAL** ——— p 38-74

Summary

Assembly guide / Patient lifts : NAUSIFLY 4

• Presentation	40
• Contraindications	40
• Contents	40
• Assembly instructions	41-43
• General safety instructions	44
• Safety instructions before first use	45-47
• Using the device	48-49
• Using the motor	50-51
• Using the sling	52-56
• Storage	57-58
• Reuse	59
• Labeling	59-61
• Cleaning and maintenance	62
• Preventive maintenance and safety checks	63-65
• Technical and dimensional specifications	66-67
• Spare parts	68
• Motor technical specifications	69-70
• Troubleshooting guide	71
• Warranty	72

NAUSIFLY 4

NAUSICAA
MEDICAL

NAUSIFLY 4 patient lift is a Class I medical device in accordance with Regulation (EU) 2017/745. ALL OUR DEVICES COMPLY WITH THE NF EN ISO 10535: 2021 STANDARD.

Presentation

NAUSIFLY 4 is a patient lift designed for use in hospitals, clinics, retirement homes, nursing homes, and at home, and must be operated by a caregiver.

It is intended for the transfer and transport of patients with reduced mobility due to illness or disability and must be used with an appropriate sling.

NAUSIFLY 4 Standard is designed for use with patients weighing up to 200 kg. NAUSIFLY 4 Compact is designed for use with patients weighing up to 160 kg.

The climatic conditions for using NAUSIFLY 4 are as follows: ambient temperature from 0°C to 40°C, humidity from 20% to 80%, and air pressure from 700 to 1060 hPa. The device can be used in the bathroom (without contact with water) or toilet. Patient lifts allow the patient to be transferred passively from a bed, recliner, toilet, or wheelchair.

Contraindications

Diseases such as osteogenesis imperfecta, osteoporosis, spinal damage, mental abnormalities, or epileptic seizures can be contraindications.

Contents

The patient lift has already been inspected to ensure that it is free of defects and that nothing is missing. Nevertheless, please check the product immediately upon receipt for any damage that may have occurred during transport. Use the delivery note to check that all items are present and that the delivery is complete.

Box	Contents	Quantity
NAUSIFLY 4 Compact and Standard	NAUSIFLY 4 Compact and Standard	1
	Control block	1
	Removable battery	1
	Remote control	1
	QR code (user manual)	1

Assembly instructions


BEFORE USING YOUR PATIENT LIFT, YOU MUST CHECK THAT:

- The lift's feet open and close correctly.
- The wheels turn and roll normally.
- The rear wheel brakes are working correctly.
- The sling bar rotates and swings correctly.
- There is no wear or deformation on the sling bar hooks.

CAUTION:


- Assembly is carried out with the brakes locked and the feet tightened.
- Do not trap the remote control wire between two metal parts.

Remove the device from its box.

Step 1:

1. Slightly loosen the front knob
2. Remove the rear knob

Step 2:


1. Raise the mast

Assembly instructions

Step 3:

1. Insert the rear knob into the lower part and tighten it firmly.
2. Tighten the top knob firmly

Step 4:

1. Loosen the knob behind the flail and the notch outlet
2. Once the flail is released, tighten the knob on the flail so that it does not get lost.

The device is ready for use.

Assembly instructions

To fold the device, follow these steps in reverse order.

USAGE TIPS:

- Your patient lift is designed to transfer patients and should not be used for any other purpose.
- Check that the patient's weight does not exceed the maximum weight that the patient lift can support.
- Operate the patient lift by pushing the handlebar, never by pushing the patient.
- Handle the lift with care when transferring a patient and at a speed appropriate to the situation.
- Move the lift on flat, smooth surfaces. It is not recommended for use on slopes greater than 5°: if you must move on a ramp, it is advisable to have a second person assist you.
- It can be used in damp rooms such as bathrooms and toilets. It is not designed for use in the shower.
- Never recharge the batteries of a patient lift near a bath or shower.

General safety instructions

1. Use the patient lift only for its intended purpose, in compliance with current legislation on medical devices, legal requirements relating to them, and regulations on occupational safety and accident prevention.
2. Please note that the patient lift is a medical device, so its user is required to comply with the directive on the use of medical devices.
3. The electrical installation requirements for the room or area where the patient lift is used must comply with the current state of the art
4. Only use the patient lift after receiving training in its use and with full knowledge of how it works.
5. Before use, read the entire user manual to avoid damage due to improper handling or exposure to hazards. The user manual contains important information and notes necessary for using the patient lift.
6. Use the patient lift only in accordance with these instructions for use. They can be viewed at any time using the QR code on the mast. Keep the instructions for use in a safe place for future reference in case of questions. Attach these instructions for use to the patient lift in case of a change of ownership.
7. Before use, check that the patient lift and its accessories are in good working order and in perfect condition.
8. Before using the lift with other medical or non-medical devices, check that the combination of these products is permitted and that they can be used together safely.
9. Assembly, commissioning, maintenance, and repair of the patient lift must only be carried out by competent specialists.
10. It is the responsibility of the user or operator to ensure (through appropriate measures and instructions) that mechanical stress on the charging cable (through bending, pulling, shearing, crushing) is prevented during charging or cleaning of the device. This also applies to the electrical cables of other devices used with the patient lift.
11. Observe the start-up time and maximum permissible load. These values must not be exceeded, otherwise safe operation can no longer be guaranteed.
12. Do not expose the patient lift to direct sunlight, heat, or humidity.
13. Ensure that no moisture enters the electrical system.
14. Avoid mechanical stress on the electrical cords used. Pulling, bending, or crushing the electrical cords may damage them.
15. Recharge batteries in a well-ventilated area.
16. Electromagnetic or other interference between the patient lift and other devices cannot be ruled out. If there is a risk of such interference, the source of the interference must be moved away or the patient lift must not be used.
17. Interference caused by the use of portable communication devices cannot be completely ruled out. For this reason, a safety distance of at least 3.3 m must be maintained to ensure the safe operation of the patient lift.
18. Do not allow children to stand unattended near the lift.
19. The lift must no longer be used if abnormal noises, damage, or any other malfunction occurs. In this case, do not connect the lift to the charger, but inform NAUSICAA MEDICAL;
20. If damaged or defective, the patient lift must no longer be used and must not be connected to the mains. Inform the dealer so that they can remedy the defect or fault.
21. Any serious incident occurring in connection with the device should be reported to NAUSICAA MEDICAL and to the competent authority of the Member State in which the user and/or patient is established.

Safety instructions before first use


Follow these safety instructions before each patient transfer:

- The caregiver must have the necessary knowledge to select and use an appropriate lifting strap.
- Check that the lift beam is compatible with the strap before using the lift.
- Check the maximum load capacity of the device and the sling.
- Check the size and model of the sling according to the patient's body type.
- Check the condition of the sling before each use. The sling must not have any tears in the fabric or damaged seams.
Check that the correct hooks have been attached. All sling buckles have 3 different levels: long (green) - medium (yellow) - short (red). Each pair of buckles must only have the following hook combination: long (green), medium (yellow), short (red).
- Check that all loops are attached to the hooks on the beam.
- Lock the wheels of the wheelchair or medical bed so that the patient can be lifted and lowered safely. Leave the casters of the patient lift in the unbraked position.
- Ensure that the transfer distance is as short as possible and never leave the patient suspended unattended.
- Observe the patient's behavior during the transfer. Sudden movements by the patient or obstacles can cause danger.
- Lift the patient to the required height.
- Keep the lift slings away from intense heat or open flames; they are not fireproof.
- Before each use, it is important to check the condition of the entire transfer device (lift + strap): seams, condition of straps and fabric.
- It is strictly forbidden to move and transport the patient while suspended when a weighing scale is mounted on the beam.

Safety instructions before first use

Before use, check that the remote control and the actuator are correctly connected to the control unit (page 32).
Check the condition of the control box and the battery.

Check that the cylinder is properly secured.

Check that the battery is loaded using the indicator on the battery. Press and hold the BATT button to display the charge level.

Safety instructions before first use

Braking system:

The brakes are an important safety feature of the lift. The rear wheels are equipped with brakes that allow you to lock the wheel by pressing your foot on the red pedal on the wheel until it stops (wheel with brake applied). To unlock the wheels, push the red pedal upward (wheel with brake released).

When locking your lift, always ensure that both wheels are locked.

Be careful when using the device on a sloping surface; both rear wheels must be locked to prevent the device from moving or tipping over.

Front wheel without brakes

Back wheel with brakes

Using the device

Patient lifts enable patients to be transferred from a bed, chair, toilet, wheelchair, or floor in a simple and efficient manner.

In healthcare premises, they offer caregivers the most appropriate solution for basic patient transfer tasks: lifting and repositioning.

TRANSFER FROM A CHAIR

non-contractual picture

TRANSFER FROM A BED

non-contractual picture

Patient lifts are designed for transferring and transporting patients with reduced mobility due to illness or disability. They must be used by a caregiver with an appropriate sling.

It is important to select a sling that is suitable for the patient based on their condition and body type. Using the right equipment significantly reduces the risk of injury to both the patient and the caregiver.

The use of a patient lift requires an assessment of the patient's physical abilities.

The following transfers can be performed with Nausicaa Medical patient lifts:


- Transferring from a chair to a bed
- Transferring from a bed to a chair
- Transferring the patient from the floor

Using the device


1. If the patient is lying on their back, turn them onto their side so that their back is facing you.
2. Fold the strap in half lengthwise.
3. Place the strap with the folded side on the patient's back. When doing so, the NAUSICAA logo and labels should be facing down. Make sure that the lower edge of the rear part of the strap is on the coccyx and that the upper edge is on the patient's shoulders.
4. Turn the patient onto the strap on the other side.
5. Pull the folded half of the strap under the patient and position it correctly.
6. Turn the patient back onto their back.
7. The patient is lying correctly on the strap if their back is completely on the back part of the strap and the leg straps are next to their thighs.
8. Now lift the backrest of the medical bed until the patient is almost sitting upright.
9. Fold the two leg supports from the outside inwards around the patient's thighs.
10. Position your patient lift so that the sling hooks are slightly above the patient's eye level. Make sure there is sufficient space between the patient's head and the patient lift sling.
11. Before attaching them to the spreader bar, make sure that the shoulder hooks and leg hooks are at the same height and symmetrical.
12. Make sure that the shoulder hooks of the strap are correctly positioned at the outer hooks of the spreader bar.
13. Then, cross the leg hooks over each other and attach them to the inner hooks.
14. Raise the lift arm of the patient lift until the shoulder and leg hooks are taut. Now check that the sling is attached correctly.
15. You can now lift the patient. Use the control handle located on the upper edge of the rear part of the sling for easier positioning.
16. Perform the reverse operation to remove the support element.

Use of motorization

Remote control for mechanical ECP version

Remote control for electrical ECP version

Use of motorization

Turning the control box on/off:

- To turn on the control box, turn the emergency stop button a quarter turn to the right. The LCD screen will then display the battery charge level.
- Press the emergency stop button after use or in case of imminent danger.

Raising/lowering:

- Press and hold the raise or lower button to activate the arm. Releasing the button stops the movement and stops the arm when it reaches the upper or lower limit.

Opening/closing the legs (ECP version):

- Pressing and holding the open or close button activates the legs of the device. Releasing the button stops the movement and prevents the device from reaching the open or closed end position.

Charging the battery:

- The battery pack is charged directly via the control box. To do this, plug the power cord (a) into the power outlet. Wall charger optional.

Removing/replacing the battery:

Wall chargers are sold as optional equipment.

To remove the battery, press the gray button and pull the unit upward to remove it.

To replace the battery, hold down the gray button and position the unit starting from the bottom, then release the gray button to snap the top part into place. release the gray button to snap the top into place.

Using the sling

Before using a patient lift, the following parameters must be taken into account depending on the situation: the patient's physical disability, medical condition, and general body type.

When selecting a sling, the caregiver must consider the following three factors:

- the patient's height
- the patient's weight
- the patient's waist circumference for slings with chest straps

All our models are Class I medical devices. They are adaptable and compatible with all patient lifts equipped with an equivalent attachment system (certificate of compatibility available on request).

FABRICS

Depending on the sling model, we offer two types of polyester fabric: JERSEY : very easy to handle, provides optimal support and comfort.

Mesh : ideal for bathing and washing, prevents maceration and is easy to dry.

Poids / Taille	1m50	1m55	1m60	1m65	1m70	1m75	1m80	1m85	1m90 + 1m90
40 kg	S	S	S	S	S	S	M	M	M
45 kg	S	S	S	S	S	S	M	M	M
50 kg	S	S	S	S	S	S	M	M	M
55 kg	S	S	S	S	S	S	M	M	M
60 kg	M	M	M	S	S	S	M	M	M
65 kg	M	M	M	S	S	S	M	M	M
70 kg	M	M	M	M	M	M	M	M	M
75 kg	L	L	M	M	M	M	M	M	M
80 kg	L	L	L	L	M	M	M	M	M
85 kg	L	L	L	L	M	M	M	M	M
90 kg	L	L	L	L	L	L	M	M	M
95 kg	XL	XL	XL	XL	L	L	M	M	M
100 kg	XL	XL	XL	XL	L	L	L	L	L
110 kg	XL	XL	XL	XL	L	L	L	L	L
120 kg	XL	XL	XL	XL	XL	XL	L	L	L
130 kg	XL	XL	XL	XL	XL	XL	XL	XL	XL
160 kg	2 XL	2 XL	2 XL	2 XL	2 XL	2 XL	XL	XL	XL
190 kg	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL	XL	XL	XL
220 kg	3 XL	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL	2 XL	XL
250 kg	4 XL	4 XL	3 XL	3 XL	3 XL	3 XL	3 XL	2 XL	2 XL
280 kg	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	3 XL	3 XL	3 XL
320 kg	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL	4 XL

IMPORTANT RECOMMENDATIONS

- To ensure maximum effectiveness when using products in this range, it is essential to:
 - choose the appropriate size for the patient
 - adjust the product to fit the patient as well as possible
- These products must not come into direct contact with broken skin.
- Read the instructions for use before using the slings and carry out the checks mentioned before each use.
- Before transferring a patient with the standing frame, ensure that the sling is securely attached to the four hooks on the lift beam.
- Sling lifespan: 3 years.

Using the sling

Read the user manual before using the slings.

The backrest straps must be hooked onto the same color.

Never adjust the sling hook settings on the lift beam when a person is seated in the sling.

Before transferring the patient, ensure that the strap is securely attached to the beam.

ISO 10535: 2021 STANDARD ON PERIODIC INSPECTION OF STRAPS

B.2.4 Inspection of body support elements (= slings)

The flexible support element should be checked regularly at the frequency specified by the manufacturer, at least twice a year.

More frequent checks may be necessary when a flexible support element is used or cleaned more frequently than normal.

Checks should be carried out by a person with the necessary qualifications and familiar with the model of the support element, as well as how it is used and maintained.

The inspection should detect signs of deterioration, wear, or potential failure, and verify the legibility of labels.

The inspection report should be kept in a secure location for reference in the event of an incident..

This report should include the following information:

- the date of the inspection;
- the identification details of the support element and its serial number;
- information concerning the condition of the support element;
- the next scheduled inspection date;
- the name, contact details, and signature of the inspector.

POSITIONING OF ATTACHMENTS

The backrest and thigh support straps must be attached to the same color.

Sitting position :

- Short high fastenings at shoulder level
- Long low fastenings, crossed between the legs

Sitting position :

- High midriff fastenings at shoulder level
- Long low fastenings

Using the sling

Semi-sitting position: :

- High fastenings at shoulder level
- Long low fastenings, crossed between the legs

For transfers in a lying position or picking up from the floor, it is necessary to always use a strap with a headrest or to add a removable headrest to the sling.

For patients without head support, it is necessary to always use a strap with a headrest or to add a removable headrest to the sling..

THIS MEDICAL DEVICE MUST BE USED OR TRAINED ON BY A HEALTHCARE PROFESSIONAL.

Lifting patients from the floor

Our Nausifly patient lift allows you to lift patients from the floor. The slings to be used are:

- Hammock slings
- U-shaped slings

When lifting a patient from the floor, you should:

- Secure the patient's neck using a rigid headrest or a headrest with fasteners.
- If you do not have a headrest, support the patient's neck with your hand.
- The patient's head should be opposite the mast to avoid impact during lifting.
- To facilitate the passage of the lift's feet, they should be spread apart, the first along the body and the second under the bend of both knees.
- **Check that the strap does not catch on the base of the cylinder at the start of the lift.**

DEMONSTRATION VIDEO

SCANNEZ
MOI

Using the sling

Patient Position

When choosing a sling, it is important to consider the different possible transfer positions in order to adapt to the patient's condition.

There are three possible positions:

- Seated position (90°)
- Semi-seated position
- Semi-reclined position

	 <p>SEATED POSITION (90°)</p>	 <p>SEMI-SEATED POSITION</p>	 <p>SEMI-RECLINED POSITION</p>
Patient Type	Patient with low dependency and adequate muscle tone	More dependent patient with a tendency to tip forward	Highly dependent patient with very low muscle tone
Scope of Use	Ideal for transferring from a chair with a fixed backrest at 90° 	Ideal for transferring from a recliner chair 	Ideal for transfers from a reclining shell chair or bed 
Recommended Slings	Eco U-shaped sling U-shaped sling Quick-adjust sling Quick transfer sling Toilet sling Ergonomic sling Wrap-around U-shaped sling Backrest sling	Eco toilet sling Mesh hammock sling Hammock sling without crossover	Hammock Sling

Using the sling

	RECOMMENDED SLINGS FOR TRANSFER		Slings recommended for washing and toilet use. Cannot remain under the patient.
	HAMMOCK SLING (Can remain under the patient)	U-SHAPED STRAP (Cannot remain under the patient)	
Loss of muscle tone in the lower limbs with maintenance of the torso and head	Hammock sling	U-shaped slings Quick-adjust U-shaped slings	Toilet sling Comfort toilet sling Quick transfer U-shaped sling
Loss of muscle tone in the lower limbs without maintaining the torso and with maintaining the head	Hammock sling	U-shaped slings Quick-adjust U-shaped slings	Toilet sling Comfort toilet sling
Loss of muscle tone in the lower limbs without support for the torso or head	Hammock straps with headrest option	U-shaped slings with backrest Quick-adjust U-shaped slings with headrest option Ergonomic U-shaped sling	
Amputated patient	Hammock strap without crossover		
Restless/anxious patient	Hammock strap without crossover	Wraparound U-shaped slings U-shaped sling with backrest	
Spastic patient	Hammock strap without crossover with optional headrest	U-shaped slings with backrest Quick-adjust U-shaped slings with headrest option Ergonomic U-shaped sling	Toilet sling with headrest option Comfort toilet sling Quick-adjust U-shaped sling
Physically and mentally dependent patient	Hammock sling	Wraparound U-shaped sling U-shaped slings with backrest Quick-adjust U-shaped slings with headrest option Ergonomic U-shaped sling	
Contraindication of hip constrictions	Hammock strap without crossover	Quick transfer sling	Quick transfer sling




Storage

Folding procedure:

Step 1: Lower the mast and crossbar to the cylinder stop using the remote control.


Step 2: Slightly loosen the knob on the beam.

Step 3: Pass the knob through the groove at the base of the jack support (a) and tighten it firmly, making sure that the beam remains horizontal (b).

Step 4: Slightly loosen the upper knob

Step 5: Loosen (a) and completely remove the lower knob while holding the device's mast (b).


Storage


Step 6: Lower the upper part toward the front.

Step 7: Insert the knob into the upper hole to lock the mast in a horizontal position.

Step 8: Tighten the knobs.

Step 9: Lift the device by the handle.

Device folded upright.

To unfold, follow the steps in reverse order (see pages 41 and 42).

Reuse


This product can be reused.

The lift and sling should be cleaned, disinfected, and maintained after each use. Before reuse, refer to the following chapters:

- Preventive maintenance and safety checks
- Cleaning and maintenance

Labelling

Labelling on the box

Labelling

1

 **SAFETY INSTRUCTIONS** - Before each use **NAUSICAA**
MEDICAL

- Check the perfect condition of the hoist sling.
- The sling must never be used to move the hoist.
- Check the perfect condition of the traction sling (Easylev range).
- Verify the presence and conditions of safety devices : all axis and their constrained fixations (pins, clips and bolts) as well as the tightness of the screws (see user manuals).

V01.03/23

DO NOT USE THE DEVICE IN CASE OF ANY LACK OR BAD OVERALL STATE

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CE

4

 **BATTERY MAINTENANCE** **NAUSICAA**
MEDICAL

A POORLY MAINTAINED BATTERY DETERIORATES

- Regularly charge the battery at the maximum.
- Batteries must be used and stored at temperatures between 5°C and 40°C.
- The charger should be checked if the battery stops working.
- Making small incomplete charges shortens the battery life.

<p>LEAD BATTERY</p> <ul style="list-style-type: none"> • Full charge twice a week when used at home. • Battery lifetime depends on its correct loading. 	<p>LITHIUM BATTERY</p> <ul style="list-style-type: none"> • Can be recharged as often as desired. • If the charge level indicator remains off, recharge the battery.
--	---

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Cleaning and maintenance

TARGET:

- Restore the patient lift and prevent the transmission of germs from one patient to another.
- Remove all organic contamination by mechanical or chemical means (disinfection).

BEFORE ANY MANIPULATION:

- Disconnect the power cable from the mains.
- Check that all electrical components are connected to each other.

CLEANING PROCEDURE:

- Clean surfaces with a damp cloth using a suitable detergent.
- Regular cleaning is recommended and should therefore be added to internal tasks.

DAILY MAINTENANCE / BETWEEN EACH PATIENT:

- Disconnect the power cord from the mains.
- Clean using a surface cleaning product.
- Clean the cylinder housing and cables if they have been splashed with bodily fluids, particularly urine.
- Clean and then disinfect all parts that come into direct contact with the patient.
- Specific maintenance by service providers after removal of the patient lift from the facility:
 - Bio-cleaning operation.
 - Steam clean the various flat surfaces. Change the cleaning surfaces regularly to avoid any water build-up. Steam clean parts that are difficult to access. For tubes, use steam with a microfiber cloth. Do not direct the steam directly onto electrical boxes.

COMPLETE CLEANING AND DISINFECTION:

- Disconnect the power cable from the mains.
- Clean the entire chassis of the device with a pre-disinfectant detergent.
- Clean the motor with a cloth slightly dampened with a pre-disinfectant detergent solution (Never clean the device, especially the electrical system, with a high-pressure cleaner, water hose, or similar).
- Disinfect the device (by spraying a surface disinfectant or misting a disinfectant solution, for example), making sure to follow the contact times for the product used.

MAINTENANCE AFTER CLEANING AND DISINFECTION:

- Disconnect the power cable from the mains.
- Check that the device has no visible damage.
- Check that no parts are missing.
- Check that the wheels are working properly and that there is no material obstructing their movement (hair, etc.).
- Check that the controls are working properly and that the remote control and actuator are connected to the battery pack.
- Clean the sockets and control buttons with a dry cloth, or a damp cloth if necessary (never clean the device, especially the electrical system, with a high-pressure cleaner, water hose, or similar).
- Check the integrity of the electrical cables (actuator and remote control).

CAUTION:

- **Detergents used must be pH neutral.**
- **Avoid abrasive products and solvents, as they could damage the surface of the appliance.**

Preventive maintenance and safety checks

Type of patient lift : NAUSIFLY 4

Chassis serial number :

Maintenance must be carried out by trained personnel.

Safety checks (frequency: depending on use, at least once a year):

This check is carried out visually; any signs of damage must result in replacement..

COMPLIANT NON-
COMPLIANT MAINTENANCE DATE

		COMPLIANT	NON- COMPLIANT	MAINTENANCE DATE
1	Fixing the flail (grooved shaft), condition of the plastic washer (and ball joint for the standard model)			
2	Lifting arm/mast attachment (bolted)			
3	Lower cylinder mounting (bolted) top and bottom			
4	Mast attachment no. 1 / base (two knobs)			
5	Attachment of the foot spreader pedal or electric cylinder (bolted)			
6	Wheel attachment (bolted)			
7	Fastening between the base and the legs (bolted)			
8	Condition of the structure (welds, joints, oxidation)			
9	Checking the presence and legibility of labels			
10	Attachment between the beam and the PESON (if PESON is mounted on the device).			

Checking electrical functions (frequency: depending on use, at least once a year):

This check is performed with the lift in use..


COMPLIANT NON-
COMPLIANT MAINTENANCE DATE

		COMPLIANT	NON- COMPLIANT	MAINTENANCE DATE
A	How the lifting/spreading cylinder works			
B	How the emergency stop works			
C	Battery status (charge retention)			
D	Remote control status			

Checks carried out on:	Checks carried out by:	Scheduled date of next inspection:

Preventive maintenance and safety checks

Checkpoints common to all 4 models


Preventive maintenance and safety checks

NAUSIFLY 4 LUBRICATION

- 1- Assembling the lifting arm/sledgehammer
- 2- Assembling the upper cylinder cap
- 3- Assembly of lifting arm/mast
- 4- Assembly of lower cylinder cap
- 5- Assembly of mast/base
- 6- Assembly of foot spreader pedal or electric cylinder
- 7- Assembly of wheels
- 8- Assembly of base/feet

Periodic lubrication: use a «3 in 1» product or similar

Preventive maintenance document to download from our website and complete.

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Technical and dimensional characteristics

Technical specifications

- Steel frame
- Laser-cut and folded sheet metal, tubes, and profiles
- Oven-baked epoxy paint
- Ball bearing wheels
- Total weight of Compact version: 38 kg
- Total weight of Standard version: 42 kg
- Compact base weight: 15.5 kg bottom / 22.5 kg top
- Standard mechanical base weight: 18 kg bottom and 24 kg top
- Maximum load for Compact version: 160 kg / Standard version: 200 kg
- Product lifespan: 8 years* (excluding electrical parts)


* Subject to regular maintenance

ENGLISH

Dimensional characteristics


Compact model

Dimensions

- A.** Total length: 1028 mm
- B.** Minimum usable height: 390 mm
- C.** Maximum usable height: 1690 mm
- D.** Total height: 1298 mm
- E.** Chassis height: 113 mm
- F.** Minimum base width: 567 mm
- G.** Maximum base width: 810 mm
- H.** Turning diameter: 1069 mm
- I.** Ground clearance: 63 mm

Standard model

Dimensions

- A.** Total length: 1292 mm
- B.** Minimum usable height: 355 mm
- C.** Maximum usable height: 1770 mm
- D.** Total height: 1459 mm
- E.** Chassis height: 113 mm
- F.** Minimum base width: 660 mm
- G.** Maximum base width:
 - Mechanical spacing : 960 mm
 - Electrical spacing : 960 mm
- H.** Turning diameter: 1337 mm
- I.** Ground clearance: 63 mm

Technical and dimensional characteristics

Compact folded :
HxLxI = 1158 mm x 567 mm x 420 mm

Standard folded :
HxLxI = 1356 mm x 661 mm x 435 mm

Spare parts

**For all spare parts requests, please
contact customer service:**

Stand Up lifts & Patient lifts

Phone : +33(0) 4 66 71 71 80

Fax : +33(0)4 66 71 71 81

Email : sav@nausicaa-medical.com

Technical specifications for motorization

Power supply:

- Removable battery pack
- Digital display showing battery life
- Separate wall charger (optional)
- IP65

Class II device

Control box:

- ABS casing
- Overload protection via thermal switches
- Electronic overload protection
- Emergency stop button
- Electromagnetic compatibility: compliant with EN 60601-1-2
- Compliant with EN 60601-1
- Protection class: II
- IP54

Actuator:

- Electric safety descent
- 24V low-voltage DC motor
- Power 24 V/120 VA
- Maximum thrust: 8000 N in standard version/ 6000 N in compact version
- Limit switches via contactor
- Stroke: 330 mm (compact) / 435 mm (standard)
- Noise level: less than 60 dB at 1 meter
- IP54

Remote control:

- 2-function/4-function remote control
- IP54
- Spiral cable remote control, low voltage 5V
- Control button pressure: less than 1N



Connection diagram

Technical specifications for motorization

Control box display:

The control box has a display that shows various information to the user.

A- Overload indicator: The load to be lifted is too heavy and may damage the device and endanger users. The lift must be relieved.

B- Abnormal usage time indicator: The stand-up lift has been used for too long compared to its down-time. This is a safety feature designed to protect the integrity of the lift. The expected operating ratio is 10% usage time to 90% rest time. Allow the device to rest before attempting to use it again.

C- Battery level indicator.

D- Low battery indicator. Risk of battery damage; the battery must be recharged before further use. An audible signal alerts the user when the battery reaches 25% of its charge.

E- Indicator for maintenance or replacement of the cylinder. The device has been used for a total of 600,000 seconds, or approximately 10,000 cycles. A general overhaul of the entire device is necessary.

Troubleshooting guide

Issue	Causes	Solutions
The unloaded device does not always descend.	Our devices need to have weight applied to them in order to descend.	Press down on the lifting arm while pressing the down button on the remote control.
The moving parts of the device are hard, making the device difficult to handle.	This is due to a lack of lubrication of the moving parts.	Lubricate moving parts.
The motor does not work, but a “click” is heard in the control box when the remote control is pressed.	<ol style="list-style-type: none"> 1. The batteries are dead. 2. The actuator cable is not connected. 3. The actuator cable is damaged. 4. The actuator or control box is damaged. 	<ol style="list-style-type: none"> 1. Charge the device. 2. Connect the jack cable. 3. Replace the jack. 4. Have the electrical system checked.
The actuator does not work, but there is no “click” sound in the control box when the remote control is pressed.	<ol style="list-style-type: none"> 1. The emergency stop button has been pressed. 2. The battery is dead. 3. The remote control cable is not connected. 4. The remote control cable is damaged. 5. The entire electrical system is damaged. 	<ol style="list-style-type: none"> 1. Turn the knob to unlock it. 2. Check the condition of the charger and then charge the battery. 3. Connect the remote control. 4. Replace the remote control. 5. Have the electrical system serviced.
The actuator stops.	<ol style="list-style-type: none"> 1. The battery is dead (sound system plus “low battery” message on the screen). 2. The device has been used for too long. 3. The weight lifted by the device is too heavy. 	<ol style="list-style-type: none"> 1. Check the condition of the charger and then charge the battery. 2. Leave the device to rest for a while (to protect the service life of the cylinder). 3. Reduce the weight.
The actuator no longer moves up or down.	<ol style="list-style-type: none"> 1. The control box is broken 2. The remote control is broken 	<ol style="list-style-type: none"> 1. Have the control box service 2. Replace the control box

Warranty

- **Article 1 :** NAUSICAA Médical S.A.S. guarantees this device against any manufacturing and assembly defects in its mechanical and electrical components, solely for devices used under the conditions specified by NAUSICAA Médical S.A.S.
 - The warranty covers mechanical and electrical parts.
 - This warranty, the terms and conditions of which are set out below, is valid for 5 years, except for batteries (lithium battery: 2-year warranty, lead battery: 1-year warranty).
- **Article 2 :** The warranty entitles the purchaser to free labor and replacement of parts recognized as defective at no cost.
- **Article 3 :** The cost of shipping the device, as well as all related expenses, shall be borne by the reseller. The goods shall always travel at the reseller's risk and under their responsibility. Under warranty: the cost of return shipping after repair will be borne by NAUSICAA Médical S.A.S. Out of warranty: the cost of return shipping is borne by the reseller, whether or not they accept the repair estimate
- **Article 4 :** The warranty does not apply if the claims are the result of:
 - an accident, misuse of the device, or negligence on the part of the purchaser.
 - transport of the device without adequate protection.
 - modification or alteration not approved by NAUSICAA Médical S.A.S.
 - the impact of external agents (natural disasters, fire, impact, humidity, flooding, lightning, etc.).
 - Installation and/or use in a manner that does not comply with technical and safety standards in the event that the device is operated in a country other than the country of purchase; and/or if the power supply is not suitable for the device's operating voltage.
 - Lack of routine maintenance
- **Article 5 :** The reseller may not invoke the warranty:
 - if the serial number of the device has been removed, modified, or rendered illegible.
 - if the device under warranty has been modified without the approval of NAUSICAA Médical S.A.S.
- **Article 6 :** No equipment will be loaned during the repair of defective equipment.
- **Article 7 :** Any warranty claims must be made through the reseller. If this is not possible, the purchaser may send the equipment directly to NAUSICAA Médical S.A.S. In this case, the purchaser must include a letter with the equipment, providing the reseller's contact details and a copy of the purchase invoice
- **Article 8 :** Spare parts under warranty will only be sent after consultation with the After-Sales Service of NAUSICAA Médical S.A.S.
Please note that defective spare parts must be returned to NAUSICAA Médical S.A.'s After-Sales Service, otherwise they will be invoiced one month after the parts have been sent by NAUSICAA Médical S.A.S.
- **Article 9 :** Defective parts replaced under or outside the warranty will be guaranteed for 6 months from the date of repair or shipment of spare parts.
- **Article 10 :** No reseller may unilaterally modify the terms of this warranty

NAUSICAA

MEDICAL

Export Department
Phone : +33(0)4 66 51 50 80
Email : export@nausicaa-medical.com
www.nausicaa-medical.com

Your contact :