



www.nausicaa-medical.com

This medical device is intended for patient positioning. This is not a restraint system.
Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

LARGE ABDOMINAL BELT

FOR ALL CHAIRS



TO COMPENSATE SLACKENING
OF ABDOMINAL MUSCLES TO
POSITION IN THE BED


Maximum useful area: + 230 cm
Height of the belt: + 18 cm
Buckle: 50 mm

Its large contact surface prevents any abdominal constriction.

Available in 4 sizes S, M, L, XL

Maximum useful area: 230 cm

Waist (in cm)



The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative. This medical device must be implemented or trained by a health professional.



CCA2HP-PF

v 05.02/23

Important Recommendations :

- In order to get a maximum efficiency of these sort of product, it is necessary :
- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient
- These products must not be in direct contact with a wounded skin.



LARGE ABDOMINAL BELT

HOSPITAL RANGE

Combination of polyurethan
foam and on both faces
Jersey/PU coating
100% polyester bias
100% polyester strap
Nylon buckle

- Place the belt around the patient's waist, smooth face outwards.
- Put the pointed end through the slit.
- Pass the strap around the backrest of the chair or around the mattress.
- Adjust the length of the straps through the Nylon buckles.
- Used in bed, the degree of restraint varies depending on the setting :
 - Slit on the stomach : kept on the back without any possibility for the patient to reach the buckles.
 - Slit in the back : possibility to move laterally.