NAUSICAA

MEDICAL

www.nausicaa-medical.com

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PUMP : NAUSIFLOW 100 MAXI 2

FOR PATIENTS FROM 30 To 450 Kg*
CLASS 1 MEDICAL DEVICE

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Content

User Manual / Pump : NAUSIFLOW 100 MAXI 2

1. Introduction	29
2. Box content	29
3. Information and safety instructions	30
3.1 Symbols	30
3.2 Caution	31
3.3 Safety instructions	32
3.4 Precautions measures	32
3.5 Electro Magnetic Compatibility (EMC)	32
4. Target groups	32
4.1 Target patients	32
4.2 Users	32
5. Operation principles	33
6. Clinical principles	33
6.1 Purpose	33
6.2 Expected clinical benefit	33
6.3 Indications	33
6.4 Contraindications	34
6.5 Other measures required	34
6.6 Side effects	34
6.7 Manual control	34
7. Product introduction	35
8. Installation	36
9. Initial start	36
9.1 Precautions & restrictions on use	37
10. Settings	37
10.1 Comfort settings	37
10.2 Cycle duration	38
10.3 Modes setting	38
10.3.1 Dynamic mode (alternating 1 sur 2)	38
10.3.2 Static mode	38
10.3.3 Care mode	39
10.3.4 Patient in sitting position	39

11. Uninstallation	40
12. Alarms	40
12.1 LOW PRESSURE alarm	40
12.2 Power outage alarm	40
13. Maintenance	41
13.1 Daily or weekly maintenance	41
13.1.1 Process	41
14. Maintenance before storage	42
14.1 General information	42
14.2 Process	42
15. Storage	43
16. Required climatic conditions	43
17. Compatible mattresses	43
18. Troubleshooting	44
19. Symbols (labels and/or packaging)	45
20. Technical features	46
21. Disposable	46
22. Lifetime	46
23. Warranty	47

Introduction

1. Introduction

Nausicaa Médical thanks you for choosing this pump.

To avoid any damage or danger caused by misuse, please read this manual in its entirety before using for the first time and before any subsequent use. This manual contains important information and notes necessary for the correct use of the pump. If you have any queries about safety warnings, precautions or use, please contact your local service provider or Nausicaa Médical After Sales Service. To avoid injury or damage due to misuse, do not use the pump until all your questions have been answered.

These instructions must always accompany and/or be supplied to any person requiring the use of this pump or to the person responsible for operating it.

The life of this pump depends essentially on correct installation and correct use.

It complies with the following standards:

EN 60601-1 : 2006/A1 : 2013

EN 60601-1-2 : 2015

EN 60601-1-11 : 2015

EU Declaration of Conformity 2017/745

The term 'system' used in this manual refers to all the components, i.e.: 'this pump associated with a support'.
The term 'support' refers to the mattress topper or mattress used with this pump.

Box content

2. Box include

Please check that the contents listed below match those you received

- 1 digitally controlled pump
- 1 user manual that can be consulted using the label with the QR Code sticker on the pump.

Information and safety instructions

3. Safety information and instructions

3.1 Symbols

To make this manual easier to read and understand, here is a brief explanation of the main symbols used.

Some of these symbols and their explanations require special attention.

USER MANUAL :

It is essential to read this manual before use.

INFO :

Provides advice or information of which the user and/or patient should be aware of.

WARNING/DANGER :

Indicates an appropriate use, operation or maintenance procedure to avoid altering and/or rendering inoperative this pump or any of its components or damaging other property.

Draws attention to a potential hazard requiring appropriate procedures or practices to prevent injury.

For more information on the symbols used, see page 45.

3.2 Caution

Read and follow the safety instructions carefully. They consist of a text or a combination of a symbol and a text. The symbol used does not replace the text!

The safety instructions apply to any person using and/or intervening in any way whatsoever on this pump or one of its components.

The person who makes this pump available to the user and/or a patient must inform him/her of the potential risks and train him/her in correct use so that no additional risk results from inappropriate behaviour.

This pump must never be used unsupervised, especially when connected to the mains.

Take particular care when children, disabled people or pets are in the vicinity of the pump.

The presence of a third party is strongly recommended, particularly in the event of this pump triggering an alarm.

Use this pump only for its intended purpose. Any other use is prohibited.

Do not use with any media other than those compatible with and marketed by Nausicaa Médical.

Information and safety instructions

3.3 Safety Instructions

 This pump's electrical equipment (cable and plug) can be dangerous if handled incorrectly.

The plug must always be accessible so that this pump can be completely disconnected from the mains if necessary.

Never use the appliance if the plug and/or cable are damaged.

Do not use this pump if it has fallen on the floor or into liquid, if it is damaged, if it shows significant wear, if it makes an unusual noise, if it is not working properly or if you are unsure of its safety. Always unplug this pump immediately after use.

Repairs and/or maintenance on this pump must only be carried out by qualified personnel or by personnel who have received the necessary training. No modifications are permitted and only original parts may be used.

To guarantee the safety of this pump, the electrical components have been tested by an external, independent testing institute. However, incorrect use can still lead to the risk of a short-circuit or electric shock.

If this pump has fallen into liquid, never attempt to recover it directly without first carefully removing the plug from the mains socket.

Do not use this pump near flammable liquids or gases. Protect it from direct sunlight. Keep it away from all sources of heat (heaters, electric blankets,) and open flames (candles, etc.).

Note that this pump sucks in air from the surrounding area, so any smoke and/or excess dust may damage its internal components.

To prevent the risk of overheating, never insert or use the pump in a closed compartment (box, etc.), never cover it with anything (sheet, blanket, etc.) and never position it on a soft, unstable surface to prevent the suction holes from becoming blocked at any time.

To prevent damage to the electrical cable, avoid all moving mechanisms (bed rails, frames, etc.).

If the support is fitted with an electrical cable sheath, it must be used.

Do not use any components other than those recommended by Nausicaa Médical. The use of a non-compliant component that cannot guarantee the correct operation of this pump will automatically invalidate the legal warranty.

 It is recommended that you report to the service provider, the manufacturer and the competent authority any serious incident* that has occurred during the use of this pump or if you consider or have reason to believe that this appliance has been tampered with.

**Serious incident* means any incident which has directly or indirectly led to, or could lead to:*

(a) The death of a patient, user or other person ;

(b) the serious deterioration, temporary or permanent, of the state of health of a patient, user or other person;

Information and safety instructions

3.4 3.4 Precautionary Measures

Protect this pump from sharp objects and surfaces, including the claws or teeth of animals.

Before use or re-use, this pump must be clean and dry.
Follow the recommendations for maintenance, cleaning and disinfection (pages 41-42)

 Store this pump as indicated in this user manual (page 43).

 Observe the minimum and maximum weights supported by this pump. Failure to comply with these weights may cause serious harm to the patient. In addition, exceeding the maximum authorised weight may lead to malfunction and cause permanent damage to the pump, rendering it inoperable.

3.5 Electro Magnetic Compatibility (EMC)

This pump complies with electromagnetic compatibility standards for electro-medical equipment.

Please note that possible electromagnetic or other influences between this pump and other electrical equipment cannot be completely ruled out.

If there is a risk of mutual interference, eliminate or move away (min. 3 m) from the source of interference.

Contact NAUSICAA MEDICAL for more detailed information on EMC..

Target groups

4. Target Groups

4.1 Target Patients

The target patients are mainly adults with a height > 150 cm and a weight between 30 kg minimum and 450* kg maximum who, due to illness, injury, disability or age, are bedridden for more than 15 hours a day, have one or more pressure sores or are considered to be at risk due to a temporary or permanent deterioration in their state of health.

* depending on the medium used

4.2 Users

Users are generally healthcare professionals, who may be assisted by carers or family members.

That's why this pump is designed to be easy and convenient to use and maintain for all users. As far as use by a member of the family is concerned, it goes without saying that he or she must have received full and appropriate instructions from experienced staff or the installer and must be capable of performing the supervisory function.

Operation principles

5. Operation principles

It works by alternately inflating the cells of the support according to a specific pattern. This pattern creates a cyclical and alternating pressure movement on the different areas of the body.

Clinical principles

6. Clinical principles

This pump is one of the components of a device intended to be used for medical purposes to prevent the development of a pressure sore and/or to assist in its care.

6.1 Purpose

To reduce pressure on the areas of the body most prone to the development of a pressure sore and more specifically where there are bony prominences.

6.2 Expected clinical benefit

Maintains tissue oxygenation in anatomical areas in contact with the support by reducing the pressure applied to cutaneous and subcutaneous tissues.

6.3 Indications

This device is mainly used for the prevention and/or treatment of all types of pressure ulcers*. It is designed to reduce the incidence of pressure sores while optimising patient comfort. It is an appropriate solution for patients who are bedridden for more than 15 hours a day and get up during the day, and whose risk assessment is classified as 'medium to high*'. As a care aid, it is indicated for a pressure sore judged to be stage 4*.

All these indications are part of a clinical judgement as well as an individual protocol including repositioning, nutritional support, skin care, wound care, etc. All other aspects of care must be undertaken by a doctor as well as a senior carer.

*Depending on medical advice, Braden scale score and support used.

 Regularly check that the system is appropriate for the patient's state of health, particularly if the patient's physical condition has changed (amputation, weight gain/loss, mental changes, etc.).

Always take a holistic view of the patient's state of health.

Always take account of contraindications. If necessary, adjust the parameters of this pump or choose a more appropriate system.

Clinical principles

6.4 Contraindications

Do not use the device if the patient has one or more of the following contraindications:

- Less than 30 kg or more than 450 kg* (*depending on the support used)
- Non-stabilised fracture of the spine,
- Cervical or skeletal traction,
- Acute multiple trauma,
- Non-stabilised post-traumatic fracture

If this is the case, use is subject to the advice of a qualified doctor.

i It is also recommended not to use this system in a hyperbaric chamber or on a stretcher.

6.5 Other measures required

The system alone is not enough; other preventive measures are also essential, including:


- Changing position (at least every 2 to 3 hours) ;
- Maintaining good skin hygiene and avoiding or minimising maceration
- Regularly changing pads, especially in cases of severe incontinence;
- Observe the condition of the skin on a daily basis, or have someone do so for you;
- Check that the patient's diet is adequate and appropriate;
- Check that the patient hydrates regularly and in sufficient quantity

6.6 Side effects

There are no scientifically proven adverse effects, although there may be a slight risk that the alternating action of this pump could cause and/or promote spasticity. In such cases, the advice of a specialist doctor should be sought.

6.7 Manual control

In alternating mode, you can check that the pressure is correctly adjusted. To do this, locate the least inflated cell of the support under the patient's buttocks and try to insert a hand flat palm up (see illustration below).

You need to be able to feel slight pressure and still be able to move back and forth.

The acceptable range between the two is approximately 25 to 40 mm.

This manual checking procedure is issued by AHCP (Agency for Health Care Policy and Research).

i If you are completely unable to move back and forth, this means that the pressure is too low (healing) > **increase the pressure.**

i If you can perform this movement without any resistance, or even with extreme ease, this means that the pressure is too high > **reduce the pressure**

Product introduction

7. Product introduction

This pump can be used in short- and long-term care environments in or outside the hospital, in retirement homes, FAMs, MASs, etc. as well as at home.

Its control panel is digitally adjustable. Its electronic management system automatically maintains pressure.

Description of setting buttons

- **On/off switch (1) :**

Switches the pump on or off.

- **Comfort pressure setting (2) :**

Five possible settings to optimise patient comfort

- **Choice of Dynamic or Static mode (3) :**

Activates / Deactivates the chosen mode

- **Alarm button (4) :**

Sets audible alarm to standby

Resets alarms.


- **Set cycle time (5) :**

Choice of cycle depending on the media used and/or according to medical recommendations.

- **Care mode (6) :**

Activates/deactivates this mode.

i For safety reasons, the care mode is automatically cancelled after 30 minutes and returns to the previously selected setting.

Installation

8. Installation

Always check that the pump is not damaged before installation.

DO NOT install or use a pump that is even slightly damaged!

- Use the suspension hooks to hang the pump on the bed frame panel

i If this is impossible, suspend it on the chassis without causing vibration, BUT never place it directly on the ground.

- Place the power supply cable into the cable passage grommet of the support (if available) or position it so that does not pose any danger, that it cannot be caught, crushed, pinched, or torn.
- Plug the power supply cable plug into the wall socket..

i Make sure that the electrical specifications are in conformity with the pump

Initial start

9. Initial start

Step 1 ► Set switch (1) to the «I» position. The green indicator light turns on and a brief beep is heard. Initially, the pump enters its self-test mode upon powering on. Only after completing this analysis does it begin to supply air.

i When the pump is turned off and regardless of the reason for this shutdown, **it retains the last parameters used (comfort pressure and cycle duration)**. It restarts in static mode until the internal pressure has reached 15 mmHg of air in the used support, then it automatically switches to alternating mode.

Step 2 ► Allow the air to escape from the pump for 2 to 3 seconds, then insert the 3-way connector of the support's air supply hose into that of the pump and let the support reach the desired comfort pressure (page 37).

! **Notice that the pressurization time depends on the support used.**

i Example: it takes approximately 25 minutes for a support measuring 88x13x200 cm to become fully operational.

Initial start

i Depending on the surface you use, the Nausicaa rapid inflation pump allows you to have the device operational in ± 2 minutes.

! If the pressure of 15 mmHg is not reached within 70 minutes, a visual and audible alarm will sound to warn of a problem (page 40).

Step 3 ▶ Once the pressure has been reached, the device is ready for use. You can now place the patient on the support and adjust the desired comfort pressure.

9.1 Precautions & restrictions on use

Take all necessary precautions

- When the pump is used at home without the intervention of nursing staff or equivalent (care assistant, etc.).
- When used in the first few days following pressure sore surgery (skin graft or flap). In this case, **use the low-pressure static mode.**
- Do not use it with a support wider than 140 cm

Settings

10. Settings

10.1 Comfort settings

The pump recalling the last parameters used, it is your responsibility to adjust the comfort pressure using the  or  and as soon as the patient is lying on the support. This setting is adjusted according to their feeling, comfort level, and positioning.

i Once this setting is made, the pump will monitor pressure changes. If the pressure drops below the level set by the user, it will automatically inflate the support so that the pressure again matches the level set by the user. Once this is achieved, the pump goes into standby until the next cycle change.

i For a suitable pressure, refer to the information on page 34 (6.7 manual control) or consult the indicative table below.

10	20	30	40	50 ¹
From 30 to 40 kg	From 30 to 65 kg	From 50 to 115 kg	From 90 to 165 ¹ kg	From 150 to 450 ¹ kg

! The information provided in this table is intended to inform and guide users. It does not constitute an official recommendation. Users are advised to verify the comfort level as specified in paragraph 6, point 6.7. Nausicaa Médical disclaims any responsibility for the use or interpretation of the information presented in this table.

¹ According to the support used and/or recommended for sitting!

Settings

10.2 Cycle duration

As with comfort adjustment, the last cycle duration used is stored in memory. It will therefore be applied from the (re)start of the pump. It is up to you to set it according to the patient's condition, pathology, feeling, comfort, medical recommendations, or the dimensions of the support used.

To modify this duration, select the desired value using the 

Examples (not exhaustive):

- after pressure sore surgery > 10 minutes
- support $\geq 120 \times 13 \times 200$ cm = 20 minutes / support $\geq 120 \times 20 \times 200$ cm = 25 minutes
- night time comfort > 20 or 25 minutes

10.3 Modes setting

10.3.1 Dynamic Mode (alternating 1 over 2)

MODE DYNAMIQUE ●

The control LED lights up to indicate that this mode is active.

It is characterized by alternating pressures, thus avoiding prolonged vascular compression that could lead to tissue hypoxia.

If this one is recommended, verify that it is indeed the one selected.

Otherwise, turn it on using the button 

 **If the pressure of 15 mmHg is not reached within 70 minutes, a visual and audible alarm sounds to warn of a problem (see page 40).**

10.3.2 Static mode MODE STATIQUE ●

The control LED lights up to indicate that this mode is active.

It allows a greater distribution of the body surface on the support. This is the principle of immersion.

This immersion among other things allows for the care of people requiring using an air support but needing relative stability of the support.

This stability is often recommended in cases of immobilization, healed fracture, neurological trauma, etc.

This mode is also required to minimize pain caused by local trauma, to promote patient rest, and also during declared end-of-life situations. Depending on the patient's weight, this mode can help reduce cell crushing (clumping) when the bed backrest is > 30° inclination.

For patients who may want to return to sleeping on a traditional mattress, whether made of standard foam or memory foam, the static mode can assist during withdrawal phases before switching to the said foam =mattress.

 When using this mode  you can adjust the desired pressure using the comfort pressure adjustment buttons.  

Settings

10.3.3 Care mode MODE SOINS ●

To activate it, press the button 

Nursing care, transferring someone to or from the bed, repositioning, etc., require a firm and stable surface.

Once this mode is activated, it automatically cancels the previous one and quickly raises the support to its maximum pressure level, without alternating effects, and this for **30 minutes**. After this time, to prevent any issues, this function cancels itself and automatically reverts to the level previously set by the user.

 **For maximum efficiency, it is recommended to activate this mode BEFORE all interventions.**


 Deactivation of this mode within 30 minutes is possible. To do so, simply press the button again. This manual action also automatically restores the previous mode.

10.3.4 Patient in sitting position (angle of 60 to 90°)

When a patient needs to be positioned semi-sitting or sitting (see the diagram below), we recommend adjusting the comfort pressure to the maximum, that is **50**.

For it to be fully effective, it is essential to make this adjustment before any change in the patient's positioning.

As soon as the pump starts, you can position the patient.

Here are some recommendations for which this position is very useful:

- To assist breathing and blood circulation.
- To enable better digestion.
- To help prevent pneumonia.
- To maintain a horizontal and non-vertical visual link.
- To have a meal or swallow a drink while minimizing the risk of choking.
- To avoid pressure sores (crushing of cells in the sacral area), mainly when the bed backrest reaches an angle between 60° and 90°.
- Etc.

 Although the energetic (alternating) mode can be used with this setting, know that the static mode offers better support.

 **It is important to be aware of safety or danger status in certain zones depending on the patient's position.**

- **CERTAIN SAFETY:** patient fully lying flat (0°) up to maximum 20°.
- **PROBABLE SAFETY:** position between 20° and 30°
- **DANGER IN THE SACRAL ZONE:** position between 30° and 60°
- **DANGER IN THE SCIATIC ZONE:** position between 60° and 90°

Uninstallation

11. Uninstallation

- Turn off the pump.
- Disconnect the power cable plug from the wall socket.
- Unplug the pump support.
- Remove the pump from the bed foot panel.
- Coil the power cable (e.g., around the casing).
- Place the pump back into its original box (if available).
- Put the pump into its carrying bag.

 You can also place the pump in the center of the support and coil the whole. This helps protect the pump during transport

Alarms

12. Alarms

This pump has 2 types of alarm.
These are distinguished by their audible and visual signals.

12.1 LOW PRESSURE alarm

This is activated to warn that LOW PRESSURE has been detected (see page 20). It is distinguished by an AUDIBLE alarm emitting a beep at regular intervals and by a VISUAL alarm whose LED flashes at the same rate as the audible alarm

12.2 Power outage alarm

It comes into operation when the mains power supply is no longer in use. It has an audible alarm that beeps at regular intervals and a visual alarm whose LED flashes at the same rate as the audible alarm.

 For greater hearing comfort, the audible alarm can be deactivated by pressing the reset button. 
If the problem is not solved, only the LED continues to flash.
When the problem has been solved, reactivate the alarms by pressing the reset button again. 
After a power cut, the system restarts with the settings used before the power cut.
before the power cut.

 **If the audible alarm has been deactivated and the problem has not been resolved after 30 minutes, it will reactivate automatically.**

Maintenance

13. Maintenance

It is **MANDATORY** to clean and disinfect the pump in accordance with the applicable national and/or local directives on the cleaning and disinfection of medical devices. This section describes the procedures to be followed. It is important to follow them.

-  To maintain maximum hygiene, the pump housing should be cleaned at least once a week.
-  When daily or weekly maintenance is carried out on this pump, it is not necessary to disconnect the pump support even if the pump is switched off.

 **The pump is only protected against liquid splashes. For safety reasons, when cleaning and disinfecting, the pump must always be disconnected from the mains supply to avoid any risk of electric shock, fire or burns, and consequently to prevent malfunction.**

 **As the pump can carry germs, we recommend that you to wear protective gloves BEFORE carrying out any handling, cleaning or disinfecting treatment.**

13.1 Daily or weekly maintenance

13.1.1 Process

- Switch off the pump (mains light off), disconnect it from the mains and **leave the bracket connected to the pump.**
- Check the pump for external damage.
- Soak a wipe, clean cloth or single-use paper in a detergent/disinfectant solution and clean the outside of the casing. If necessary, leave the solution to act as recommended.

 **Take care not to apply too much pressure to the control panel. otherwise it could be damaged.**

- Make sure that all areas are clean, then wipe the housing with a clean, dry cloth and leave to dry.
- Visually check the condition of the filter. There should be no excess dust. If there is little excess (a slightly grey filter), you can remove it by removing the filter from its position and blowing it out using a dust-removal spray or low-pressure compressed air. If the filter is more heavily soiled (dark grey or black), it must be replaced.
- Once the pump has been cleaned, disinfected, dried and inspected, make sure that there is no residual moisture on the pins of the then plug it in and test it to ensure that it is working properly.

Maintenance before storage

14. Maintenance before storage

14.1 General information

It is recommended that any faulty components be replaced and that any servicing be carried out before storing the pump. This work must be carried out by qualified personnel.

14.2 Process

 Before any handling and/or intervention, wear protective gloves and carry out an initial surface disinfection.

- Check the integrity of the power supply cable. There must be no visible damage or signs of significant wear.
- Check the plug and its pins.
- Check that the casing is completely intact.
- If the pump is in order, move it to the cleaning area and place it on a stable, firm work surface.
- Soak a wipe, clean cloth or single-use paper in a detergent solution and clean the entire casing, including the power supply cable except for the plug.
- Make sure that the entire casing is clean.
- Disinfect the enclosure and power cable using the same procedure as for cleaning.

 **Be careful not to press too hard on the control panel to avoid damaging it. Strictly adhere to the usage recommendations of the disinfectant solution.**

- Once the remanence time has elapsed, if necessary, wipe off any excess solution remaining on the two components (casing and electric cable) with a clean, dry, single-use cloth and leave to dry.
- Replace the filter in any condition. To do this, remove the filter holder from the back of the pump, remove the old filter, clean and dry the filter holder, then insert the new filter and replace the filter holder.
- Once the pump has been cleaned, disinfected and dried, make sure that no residual moisture remains.
- Check that there are no defects in the internal components and that the fastening system is intact. If necessary, repair and/or replace defective components.
- Connect and test the pump to ensure it is working properly.
- Once operation has been checked, unplug the pump, put it back in its box (if available) and affix a visible identification label to ensure that it has been cleaned, disinfected and is working properly, then put it away.

 **DURING THESE INSPECTION PHASES, ENSURE THAT:**
DO NOT immerse or soak the pump.
DO NOT spray liquid directly onto the surface of the pump to prevent any damage caused by liquid penetrating inside the casing.
DO NOT use cleaning solutions based on hypochlorite, phenolic compounds, corrosive products, solvents such as acetone, ether, or other similar products such as iodine alcohol, etc., as this could damage the casing.
DO NOT use abrasive materials.
DO NOT leave excess cleaning and/or decontamination solution on the control panel, as this could cause damage if it penetrates.
DO NOT place the pump directly on the ground.
DO NOT wash the filter

Storage

15. Storage

We recommend :

- Store the pump in a transport bag or in the original packaging.
- Handle the packaged pump with care.
- Do not place heavy objects on it.
- Protect it from direct sunlight.
- Respect the climatic conditions for storage

Required climatic conditions

16. Required climatic conditions

Conditions	Temperatures (°C)	Humidity level	Athmospheric Pressure (hPa)
When in use	From 10°C to 40°C	From 30% to 75% (without condensation)	From 700 hPa to 1.060 hPa
Storage/shipping	From 0°C to 40°C		

Compatible mattresses

17. Compatible mattresses

Height in cm	Dimensions in cm (L x la)	Foam Layer	Max. User Weight
13	70 x 200	YES	150 kg
	88 x 200		
	120 x 200		
20 (full air)	85 x 200	NO	200 kg
	120 x 200		
25 (full air)	120 x 200	NO	450 kg
	140 x 200		

Troubleshooting

18. Trouble shooting

The following table lists some problems, possible causes, and solutions. If the solutions in the table do not resolve the issue, contact your distributor or Nausicaa Medical after-sales service.

Failure and/or alarm	Possible cause	Solution
Pump is not fonctionning	The plug is not inserted into the wall socket or only partially inserted	Fully insert the plug into the wall socket
	The power outlet is faulty or the power supply is off	Change the outlet or restore power
	The on/off button is not pressed	Press the on/off button
	The power cable is damaged	Replace the power cable
	1 or 2 fuses are blown	Replace fuse(s)
	If the power indicator is on and none of the causes listed above are observed, contact your local provider or Nausicaa Medical after-sales service.	
Pump is not fonctionning	The static mode is selected	Select dynamic mode
	Internal problem (mechanical or other)	Contact your local provider or Nausicaa Médical after-sales service
High permanent pressure	Internal problem (mechanical or other)	Contact your local provider or Nausicaa Médical after-sales service
No or very little air from the pump	Dirty air filter Faulty pump	Check the air filter Replace pump
Low pressure in the supporting cells	The set pressure is too low considering the patient's weight and/or position (sitting). Internal problem (mechanical or other)	Increase the pressure Adjust according to the patient's position Contact your local provider or Nausicaa Médical after-sales service
Low Pressure Alarm	CPR open Air leak at the connecto Air leak in cell Internal problem (mechanical or other)	Close the CPR Check the connector Check cells and connections Contact your local provider or Nausicaa Médical after-sales service

Symbols (labels and/or packaging)

19. Symbols (labels and/or packaging)

	Indicate that the user manual must be read		Indicate an important warning
	Advice or information that the user and/or the patient should be aware of		Indicates that the user manual must be read before using the product
	Indicate that the product is a medical device compliant with the EU Medical Devices Regulation 2017/745		Indicates that the product complies with European Union legislation
	Unique device identification code		Manufacturer contact details
	Manufacturing date		Serial number
	Article reference		Complies with IEC 60601-1 type BF standard
IP21	Electrical Protection index		Indicates the type of fuse
	Electrical insulation class II: double insulation		Do not dispose with household waste
	Do not use chlorinated bleaching agents (Cl)		User's weight

Technical features

20. Technical features

Model	NAUSIFLOW 100 MAXI 2
Reference	NA100MAXI2-COMP-PM
MDR Class	Classe I
Basic UDI Number	37014294NAUSIFLOWAUTOMGE
UDI-DI Number	03701429409142
Electrical Protection Class	 Classe II, Type BF, (No AP/APG)
Noise Nuisance	27-28 dB (maximum flow rate at 1 m) (ISO 7779)
Mains Voltage	220 V
Mains Frequency	50 Hz
Fuse Rating	250 V, T1A
Protection Index	IP21
Air Flow	14 L/min (double pump).
Pressure Range	15 à 50 mmHg
Cycle Duration	10, 15 & 20 minutes
Outside pump material	ABS plastic
Dimensions	(L) 40 mm x (L) 14 mm x (H) 20 mm
Weight	3.85 kg

Disposable

21. Disposable

Respect national regulations for the removal of one of the components of this pump or for its complete disposal.

Lifetime

22. Lifetime

The lifetime indicated below depends on the frequency and intensity of use, cleaning, disinfection, transport, compliance with annual maintenance, etc.

Product	Lifetime
Pump NAUSIFLOW 100 MAXI 2	+ five (5) years

 Indicating a specific lifetime does not mean an additional warranty

Warranty

23. Warranty

Article 1: NAUSICAA Médical S.A. guarantees this appliance against all manufacturing and assembly defects in its components, and only for appliances used under the conditions stipulated by NAUSICAA Médical S.A. This guarantee, the terms of which are defined below, is valid for 24 months from the date of first purchase from NAUSICAA Médical S.A.

Article 2: The warranty entitles the purchaser to free labour and the replacement, free of charge, of any parts found to be defective.

Article 3 : The reseller is responsible for the carriage of the appliance and all related costs. Goods always travel at the risk and under the responsibility of the retailer. Under warranty: NAUSICAA Médical S.A. will be responsible for the cost of returning the goods after servicing. Out of guarantee: the return costs are at the retailer's expense, whether or not he accepts the repair estimate.

Article 4: The guarantee does not apply if the claims are the result of :

- An accident, misuse of the appliance or negligence on the part of the purchaser.
- Transport of the appliance without adequate protection.
- Modifications or alterations not approved by NAUSICAA Médical S.A.
- The impact of external agents (natural disasters, fire, impact, humidity, flooding, lightning, etc.).
- Installation and/or use that does not comply with technical and safety standards, if the appliance is to be used in a country other than the country of purchase; and/or if the electrical supply is not suitable for the voltage at which the appliance is to be used.
- Failure to carry out routine maintenance.

Article 5: The retailer may not invoke the benefit of the guarantee:

- If the serial number of the appliance has been removed, modified or rendered illegible.
- If the equipment under guarantee has been modified without the approval of NAUSICAA Médical S.A.

Article 6: No equipment will be loaned during the repair of defective equipment.

Article 7: All warranty claims must be made through the retailer. If this is not possible, the purchaser may send the equipment directly to NAUSICAA Médical S.A. In this case, the purchaser must indicate the retailer's contact details and a copy of the purchase invoice in a letter enclosed with the equipment.

Article 8: Spare parts will only be sent under guarantee after consultation with NAUSICAA Médical S.A.'s After Sales Service.

It should be noted that defective spare parts must be returned to the NAUSICAA Médical S.A. After-Sales Service within 15 days, failing which they will be invoiced 1 month after the parts have been sent by NAUSICAA Médical S.A.

Article 9: Defective parts changed under or outside the guarantee are guaranteed for 6 months from the date of repair or dispatch of spare parts.

Article 10: No retailer may unilaterally modify the terms of this warranty

NAUSICAA

MEDICAL

Customer Service

Phone : 04 66 51 50 80

Fax : 04 66 51 50 47

Mail : contact@nausicaa-medical.com

www.nausicaa-medical.com

Dealer:

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© Nausicaa Médical S.A.