# NAUSICAA

## MEDICAL

[www.nausicaa-medical.com](http://www.nausicaa-medical.com)

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## PUMP : **NAUSIFLOW 2S**

FOR PATIENTS FROM 30 TO 200\* Kg  
CLASS I MEDICAL DEVICE

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EN

USER MANUAL  
PUMP: **NAUSIFLOW 2S** ——— p 24-48

© Nausicaa Médical S.A. Design policy and copyright:

As our policy is constantly evolving, we reserve the right to modify designs without prior notice. The contents of this publication may not be copied in whole or in part without the express permission of NAUSICAA MEDICAL.

FRANCAIS

ENGLISH

# Summary

## User Manual / Pump : NAUSIFLOW 2S

1. Introduction	28
2. Box contents	28
3. Safety informations and instructions	29
3.1 Symbols	29
3.2 Important	29
3.3 Safety instructions	30
3.4 Precautions	31
3.5 Electro Magnetic Compatibility (EMC)	31
4. Targets groups	31
4.1 Patient targets	31
4.2 Users	31
5. Operating principle	32
6. Clinical principles	32
6.1 Purpose	32
6.2 Expected clinical benefit	32
6.3 Indications	32
6.4 Contraindications	33
6.5 Other measures required	33
6.6 Side effect	33
6.7 Manual control	33
7. Presentation	34
7.1 Description of adjustment buttons	34
8. Installation	35
9. Getting started	35
9.1 Precautions and restrictions of use	36
10. Settings	36
10.1 Comfort settings	36
10.2 Cycle time	37
10.3 Mode settings	37
10.3.1 Mode dynamic	37
10.3.2 Mode static	37
10.3.3 Sitting position	38

11. Uninstall	39
12. Low pressure alarm	39
13. Maintenance	39
13.1 Daily or weekly maintenance	40
13.1.1 Process	40
14. Maintenance before storage	40
14.1 General information	40
14.2 Process	40
15. Storage	42
16. Climatic requirements	42
17. Compatible mattresses	42
18. Troubleshooting	43
19. Symbols (product label and/or packaging)	44
20. Technical features	45
21. Disposal	45
22. Lifespan	45
23. Warranty	46

# Introduction

## 1. Introduction

Dear customer,

Nausicaa Médical thanks you for choosing this pump.

In order to avoid any damage or danger caused by misuse, please read this manual in its entirety before using your pump for the first time, and before using it at any time in the future. This manual contains important information and notes necessary for the correct use of the pump. If you have any questions about safety warnings, precautions or use, please contact your local service provider or Nausicaa Médical Customer Service. To avoid injury or damage due to misuse, do not use this pump until all your questions have been answered.

This manual must always accompany and/or be supplied to any person requiring the use of this pump, or to the person responsible for operating it.

The life of this pump depends essentially on correct installation and use.

It complies with the following standard(s):

EN 60601-1 : 2006/A1 : 2013

EN 60601-1-2 : 2015

EN 60601-1-11 : 2015

EU Declaration of Conformity 2017/745

The term “support” used in this manual refers to the mattress topper or mattress used. The term “system” refers to the set of components enabling its use, i.e. “this pump and a support”.

## 2. Box contents

Please check that the contents listed below correspond to those you received.

1 manually-adjustable pump

1 user manual with QR Code label on pump.

# Safety informations and instructions

## 3. Safety informations and instructions

### 3.1 Symbols

To make this manual easier to read and understand, here is a brief explanation of the main symbols used.

Some of these symbols and their explanations require special attention.

User manual :

Please read this manual before use.

INFO :

Provides advice or information you must be aware of as user and/or patient

**CAUTION/DANGER:**

**Indicates proper use, operation or maintenance procedures to avoid altering and/or rendering inoperative this pump or any of its components, or damaging other property.**

**Warning regarding a potential issue requiring appropriate procedures or practices to prevent injury.**

### 3.2 Important

Read and follow safety instructions carefully. They consist of text or a combination of text and symbols. **The symbol used does not replace the text**

For further information on the symbols used, see page 44.

The safety instructions apply to anyone using and/or intervening in any way with this pump or any of its components.

The person who makes this pump available to the user and/or a patient must inform him/her of the potential risks and train him/her in its correct use, so that no additional risk results from inappropriate behavior.

**This pump must never be left unattended, especially when connected to the mains.**

**Strict vigilance is required when children, disabled people or pets are in the vicinity of the pump.**

**The presence of a third party is strongly recommended, especially in the event of this pump triggering an alarm.**

**Use this pump only for its intended purpose. Any other use is prohibited.**

**DO NOT use with any media other than those compatible with and marketed by Nausicaa Médical.**

# Safety informations and instructions

## 3.3 Safety instructions

 The electrical equipment (cable and plug) of this pump can be dangerous if handled incorrectly.

The plug must always be accessible so that, if necessary, the pump can be completely disconnected from the mains.

Never use the unit if the plug and/or cable are damaged.

Do not use this pump if it has fallen on the ground or into liquid, if it is damaged, if it shows significant wear, if it makes an unusual noise, if it is not working properly or if you are unsure of its safety. Always unplug this pump immediately after use.

Repairs and/or service work on this pump may only be carried out by qualified personnel or by personnel who have received the necessary training. No modifications are permitted and only original parts may be used.

To guarantee the safety of this pump, the electrical components have been tested by an external, independent testing institute. However, incorrect use can still lead to the risk of a short-circuit or electric shock.

If this pump has fallen into liquid, never retrieve it directly without first carefully removing the plug from the mains socket.

Do not use this pump near flammable liquids or gases. Protect it from direct sunlight. Keep away from sources of heat (heaters, electric blankets, ...) and open flames (candles, etc.). Please note that this pump draws in air from the surrounding area, so any smoke and/or excessive dust may damage its internal components.

To prevent the risk of overheating, never insert or use the pump in a closed compartment (box, etc.), never cover it with anything (sheet, blanket, etc.) and never position it on a soft, unstable surface to prevent the suction holes from becoming blocked at any time.

To prevent damage to the power cable, avoid all moving mechanisms (bed rails, chassis, etc.).

If the installed support has a cable duct, it must be used.

Do not use any components other than those recommended by Nausicaa Médical. The use of any non-conforming component which cannot guarantee the correct operation of this pump will automatically invalidate the legal warranty.

 It is recommended that you report any serious incident\* involving the use of this pump, or if you consider or have reason to believe that this device has been tampered with, to the service provider, the manufacturer and the relevant authority.

\*« serious incident» means any incident which has directly or indirectly led to, or could lead to:

(a) The death of a patient, user or other person;

(b) The serious deterioration, temporary or permanent, of the state of health of a patient, user or other person;

# Safety informations and instructions

---

## 3.4 Precautions

Protect this pump from sharp objects and surfaces, including the claws and teeth of animals.

Before use or re-use, this pump must be clean and dry.

Follow the recommendations for maintenance, cleaning and disinfection (pages 39-41).

 Store the pump as described in this manual (page 42).

 Observe the minimum and maximum weights supported by this pump. Failure to do so may result in serious harm to the patient. Furthermore, exceeding the maximum weight allowed may lead to malfunction and internal damage to the pump, rendering it inoperable.

## 3.5 Electro Magnetic Compatibility (EMC)

This pump complies with electromagnetic compatibility standards for electro-medical devices.

Please note that electromagnetic and other influences between this pump and other electrical equipment cannot be completely ruled out. If there is a risk of mutual interference, remove or move away (min. 3 m) from the source of interference.

Contact NAUSICAA MEDICAL for more detailed information on EMC.

# Targets groups

---

## 4. Targets groups

### 4.1 Patient targets

Target patients are mainly adults with a height > 150 cm and a weight between 30 kg minimum and 200\* kg maximum who, due to illness, injury, disability or age, are bedridden for more than 15 hours a day, have one or more pressure sores, or are considered to be at risk due to a transient or permanent deterioration in their state of health.

\*Depending on the medium used

### 4.2 Users

Users are generally healthcare professionals, who may be assisted by caregivers or family members.

That's why this pump is designed to be easy and convenient to use and maintain, whatever the user.

In case of use by a family member, it is essential that he or she has received full and appropriate instructions from the trained staff or installer, and is able to perform the supervisory function.

# Operating principle

## 5. Operating principle

The operating principle is the alternating inflation of the cells of the support used according to a specific pattern. This pattern creates a cyclic and alternating pressure movement on the different zones of the body.

## Clinical principles

### 6. Clinical principles

This pump is one of the components of a system designed to be used for medical purposes to prevent the development of a pressure sore and/or to assist in its care.

#### 6.1 Purpose

Reduce pressure on areas of the body most prone to pressure sore development.

#### 6.2 Expected clinical benefit

Maintains tissue oxygenation in anatomical areas in contact with the support by reducing the pressure applied to cutaneous and subcutaneous tissues.

#### 6.3 Indications

This **system** is mainly recommended for the prevention and/or treatment of all types of pressure sores (according to medical advice).

It is designed to reduce the incidence of pressure sores while optimizing patient comfort. It is an appropriate solution for patients confined to bed for more than 15 hours a day and getting up during the day, whose risk assessment is classified as “medium to very high\*”. As a care aid, it is indicated for a stage 4\* pressure sore.

All these indications are part of a clinical judgment as well as an individual protocol, including repositioning, nutritional support, skin care, wound care, etc. All other aspects of care must be individually assessed. All other aspects of care must be addressed by a physician as well as an experienced caregiver.

**\*Depending on medical opinion, Braden scale score assessment and media used.**

 Regularly check the **system**'s suitability for the patient's state of health, especially if the patient's physical condition has changed (amputation, weight gain/loss, mental changes, etc.).

Always assess the patient's state of health holistically.

Always take contraindications into account. If necessary, adjust the parameters of this pump or choose a more suitable system.

# Clinical principles

## 6.4 Contraindications

Do not use the device if the patient has one or more of the following contraindications:

- Less than 30 kg or more than 200 kg\* (depending on the mattress used)
- Non-stabilized spinal fracture,
- Cervical or skeletal traction,
- Acute multiple trauma,
- Non-stabilized post-traumatic fracture.

Should this be the case, use is subject to the advice of a qualified physician.

**i** It is also recommended not to use this **system** in a hyperbaric chamber or on a stretcher.

## 6.5 Other measures required

The **system** alone is not enough; other preventive measures are also essential, including :

- Changing position (at least every 2 to 3 hours);
- Maintaining good skin hygiene and avoiding or minimizing maceration;
- Regular replacement of incontinence protection, especially in cases of severe incontinence;
- Cutaneous check on a daily basis;
- Check that nutrition is adequate and appropriate;
- Check that the patient hydrates regularly and in sufficient quantity.


## 6.6 Side effect

There are no scientifically-documented adverse effects, although there may be a slight risk that the alternating action of this pump may cause and/or promote spasticity.

In such cases, specialist medical advice should be sought.

## 6.7 Manual control

In alternating mode, you can check whether the pressure is correctly adjusted. To do so, locate the least swollen cell under the patient's buttocks and try to insert a hand flat palm-up (see illustration below).

You need to be able to feel slight pressure while still being able to move back and forth.

The acceptable range between the two is about 25 to 40 mm.

This manual verification procedure is issued by the Agency for Health Care Policy and Research (AHCPR).

**i** If you are completely unable to move back and forth, the pressure is too low (healing) > **increase the pressure.**

**i** If you can perform this movement without any resistance, or even with extreme ease, this means that the pressure is too high > **reduce the pressure.**

# Presentation

## 7. Presentation


This pump can be used in short- and long-term care environments in or outside the hospital, in retirement homes, FAMs, MASS etc., as well as at home.

This pump features manual adjustment.  
It belongs to the category of “mechanical” pumps.

### 7.1 Description of adjustment buttons

- **On/off switch (1) :**  
Starts or stops pump operation.
- **Comfort pressure setting (2) :**  
Optimizes patient comfort.
- **Static or Dynamic mode (3) :**  
Activates/deactivates the chosen mode.
- **Alarm (4) :**  
Light alarm  
Resets alarms.
- **Normal or low pressure (5) :**  
Displays current pressure.

 In the event of low pressure, a visual and audible alarm is triggered.

# Installation

## 8. Installation

Always check the pump for damage before installation.

**DO NOT install or use a pump that is even slightly damaged.**

- Hang the pump on the footboard panel using the suspension hooks.

- ⓘ If this proves impossible, hang it from the chassis, ensuring that it does not cause any vibration BUT NEVER place it on the ground.
  - Position the power cable so that it does not constitute a hazard, and cannot be caught, crushed, trapped or pulled.
- ⓘ If the support you're using has a cable gland, make sure you place the electrical cable in it.
  - Insert the power cable plug into the wall socket.
- ⓘ Make sure that the electrical specifications of your mains supply meet those required for the operation of this pump.

## Getting started

### 9. Getting started

**Step 1** - Set switch (1) to "I" position. The green LED lights up and a beep is heard. This is the audible indicator of low pressure, which is normal during start-up. Press the alarm reset button (4) to stop the beep; only the LED continues to flash. This will go out once the internal pressure has been reached.

- ⓘ At startup, dynamic mode is the default mode.

# Getting started

**Step 2** - Allow the air to escape from the pump for a few seconds, then insert the 3-way connector of the bracket's air supply hose into that of the pump, position the adjustment knob (2) on the pictogram  and allow the bracket to inflate.

 **Make sure both connectors (bracket and pump) are properly connected and secured.**

 Please note that the time required to install the patient depends on the support used. For example: a mattress measuring 88x13x200 cm takes around 30 minutes to be fully operational.

 **Using the Nausicaa rapid inflation pump, you can reach the required pressure in  $\pm 2$  minutes.**

**Step 3** - Once the pressure has been reached, the normal pressure indicator lights up. The system is now ready for use. You can now place the patient on the support and adjust the desired comfort level (page 36).

## 9.1 Precautions and restrictions of use

Take all necessary precautions

- When the pump is used at home without the intervention of nursing staff or similar (care assistant, etc.).
- When the pump is used in the first few days of a post-operative pressure sore (skin graft or flap). **In this case, use the low-pressure static mode.**
- Do not use with a support wider than 100 cm.

# Settings

## 10.1 Comfort settings

Once the patient is lying on the support, you can adjust the pressure by turning the knob (2) clockwise, depending on how the patient feels, how comfortable he or she is and how he or she is positioned.

 To obtain the correct pressure, refer to the information on page 33 (6.7 manual control) or refer to the table below.

	SOFT	STANDARD	TOUGH	
From 30 to 40 kg	From 40 to 60 kg	From 50 to 100 KG	De 90 à 150 kg	Depending on position & 200 <sup>1</sup> kg max

 The information provided in this table is intended to inform and guide users. It does not constitute an official recommendation. Users are advised to check the state of comfort as specified in paragraph 6 point 6.7. Nausicaa Médical accepts no responsibility for the use or interpretation of the information contained in this table.

# Settings

 This symbol on the bottom left of the comfort pressure knob means Low Pressure, while the symbol on the knob  indicates Static Mode.

**1 Depending on the media used!**

## 10.2 Cycle time

This is set at a moderate value of 12 minutes.

## 10.3 Mode settings

### 10.3.1 Mode dynamic

It is characterized by alternating pressures, thus avoiding prolonged vascular compression that could lead to tissue hypoxia.

This is the mode generally required for most people.

As stipulated (page 35 > 9 Getting started), this is the mode that comes into operation as soon as the pump starts up. You don't need to intervene.

### 10.3.2 Mode Static

It allows the body's surface area to be more evenly distributed over the surface. This is the immersion principle.

Among other things, it enables the treatment of people who require the use of an air support, but who need relative stability of the support.

This stability is often recommended in cases of immobilization, consolidated fractures, neurological trauma, etc.

This mode is also required to minimize pain due to local trauma, to promote patient rest, and for declared end-of-life.

Depending on the patient's weight, this mode can claim to reduce cell crushing (healing) when the bed backrest inclination is  $> 30^\circ$ .

For patients who are likely to return to sleeping on a conventional mattress, whether made of standard foam or memory foam, the static mode can help to achieve weaning phases before the foam mattress is put in place.

To switch it on, press button (3) to switch from alternating to static mode; the control LED lights up to indicate that it is on. If necessary, adjust pressure by turning knob (2) clockwise. Press knob (3) again to return to previous mode.


# Settings

## 10.3.3 Sitting position

For it to be fully effective, the pictogram  must be selected **before** any changes are made to patient positioning.

To do this, turn the adjustment knob (2) clockwise.

Wait 1 to 2 minutes, then place the patient in a semi-seated or seated position, i.e. at an angle of between 30° and 90° (see photo below).

Here are just a few of the reasons why this position is so useful:

- To facilitate breathing and blood circulation.
- For better digestion.
- To help prevent pneumonia.
- To maintain a horizontal rather than vertical visual connection
- To take a meal or swallow a drink while minimizing the risk of choking.
- To restrict the talon effect (crushing of cells located in the sacral zone) as soon as the bed backrest reaches an angle > 30°.
- Etc.

 Dynamic (alternating) mode is generally used with this setting, but static mode offers a more uniform layout.

## 10.3.4 Other

For nursing care, transfer from or to bed, repositioning, etc., a firm, stable surface is recommended.

To facilitate these various actions, activate the static mode, then turn the knob (2) to the pictogram . Allow the support to inflate and present a firm surface before operating.

 **For maximum efficiency, we recommend activating this function BEFORE any intervention.**

 **To avoid any problems, don't forget to manually return to the previous mode.**

# Uninstall

---

## 11. Uninstall

- Switch off the pump.
- Remove the power cable plug from the wall socket.
- Disconnect pump bracket.
- Remove pump from footboard.
- If necessary, wind the power cable around the housing.
- Replace pump in original box (if available).
- Return the **system** to its carrying bag.

 Alternatively, place the pump in the center of the holder and roll up the wrapping. This will protect the pump during transport.

## Low pressure alarm

---

### 12. Low pressure alarm

If the pressure remains permanently low and/or the pump fails to reach or maintain the desired pressure, the AUDIBLE and VISUAL alarm is triggered.

This alarm is distinguished by a beep at regular intervals, combined with LED flashing at the same rate as the audible alarm.

For audible comfort, you can deactivate the audible alarm by pressing the reset button (4). Until the problem is solved, only the LED continues to flash. When the problem has been solved, press the reset button (4) again to reactivate the alarms.

## Maintenance

---

### 13. Maintenance

It is **MANDATORY** to clean and disinfect the pump in accordance with applicable national and/or local guidelines for cleaning and disinfecting medical devices.

This section describes the procedures to be followed. It is important to follow them.

 To maintain maximum hygiene, the pump housing should be cleaned at least once a week.

 When daily or weekly maintenance is carried out on this pump, it is not necessary to disconnect the pump support, even if the pump is switched off.

 **As the pump is only protected against liquid splashes, as a safety measure during cleaning and disinfection, it must always be disconnected from the mains to avoid any risk of electric shock, fire, burns and, consequently, to prevent malfunction.**

# Maintenance

---

 **As the pump can carry germs, we recommend that you wear protective gloves BEFORE carrying out any handling, cleaning or disinfection treatment.**

## 13.1 Daily or weekly maintenance

### 13.1.1 Process

- Switch off the pump (mains indicator off), disconnect it from the mains and leave the **bracket connected** throughout the operation.
- Check pump for external damage. Wipe, clean cloth or disposable paper in a detergent/disinfectant solution and clean the outside of the housing. If necessary, leave the solution to act as recommended.

 **Take care not to over doing as this may damage the control panel.**

- Make sure all areas are clean, then wipe the housing with a clean, dry fabric and dry.
- Visually check the condition of the filter. There should be no excess dust. If there is little excess dust (a slightly gray filter), remove the filter from its housing and clean it with a dust-removal spray or low-pressure compressed air. If the filter is more heavily soiled (dark grey or black), it must be replaced.
- Once the pump has been cleaned, disinfected, dried and inspected, make sure that no residual moisture remains on the plug pins before putting it back into service, then plug it in and test it to make sure it's working properly.

## Maintenance before storage

---

### 14. Maintenance before storage

#### 14.1 General

It is recommended to replace any defective components and to carry out complet maintenance before storing the pump. This work must be carried out by qualified personnel.

#### 14.2 Process

 Before any handling and/or intervention, wear protective gloves and carry out an initial surface disinfection.

- Check the integrity of the power supply cable. There must be no visible damage or signs of significant wear.
- Check the plug and its pins.
- Check that the casing is completely intact.
- If the pump is operational move it to the cleaning area and place it on a stable, firm working surface.

# Maintenance before storage

- Wipe, clean cloth or single-use paper in a detergent solution and clean the entire casing, including the power supply cable except for the plug.

 **Be careful not to press too hard on the control panel and electrical connections of the plug to avoid damaging it.**

- Make sure that the entire casing is clean.
- Disinfect while using the same procedure as for cleaning.

## **Strictly adhere to the disinfectant's recommendations.**

- Once the remanence time has elapsed, if necessary, wipe off any excess solution remaining on the two components (casing and electric cable) with a clean, dry, single-use cloth and leave to dry.
- Replace the filter in any condition. To do this, remove the filter holder from the back of the pump, remove the old filter, clean and dry the filter holder, then insert the new filter and replace the filter holder.
- Once the pump has been cleaned, disinfected and dried, make sure that no residual moisture remains.
- If in doubt and if the pump is out of warranty, check that there are no defects in the internal components and that the fastening system is intact. If necessary, repair and/or replace defective components.
- Connect and test the pump to ensure it is working properly.
- Once operation has been checked, unplug the pump, put it back in its box (if available) and affix a visible identification label to ensure that it has been cleaned, disinfected and is working properly, then put it away.

 **DURING THESE INSPECTION PHASES, ENSURE THAT:**

**DO NOT immerse or soak the pump.**

**DO NOT spray liquid directly onto the surface of the pump to prevent any damage caused by liquid penetrating inside the casing.**

**DO NOT use cleaning solutions based on hypochlorite, phenolic compounds, corrosive products, solvents such as acetone, ether, or other similar products such as iodine alcohol, etc., as this could damage the casing.**

**DO NOT use abrasive materials.**

**DO NOT leave excess cleaning and/or decontamination solution on the control panel, as this could cause damage if it penetrates.**

**DO NOT place the pump directly on the ground.**

**DO NOT wash the filter.**

# Storage

## 15. Storage

We recommend :

- Store the pump in a transport bag or in the original packaging.
- Handle the pump with care.
- Do not place heavy objects on it.
- Protect it from direct sunlight.
- Respect climatic storage conditions.

ENGLISH

## Climatic requirements

### 16. Climatic requirements

Conditions	Temperatures (°C)	Humidity level (%)	Atmospheric pressure (hPa)
User	From 10°C to 40°C	From 30% to 75% (non-condensing)	From 700 hPa to 1.060 hPa
Storage/ Transport	From 0°C to 40°C		

## Compatible mattresses

### 17. Compatible mattresses

Air height in cm	Dimensions in cm (L x la)	Foam base	Maximum loading weight
13	70 x 200	YES	150 kg
	88 x 200		
20 (full air)	85 x 200	NO	200 kg

# Troubleshooting

## 18. Troubleshooting

List of potential issues, their possible causes and how to solve them. If the solution below does not resolve the problem, contact your distributor or Nausicaa Médical's after-sales service.

Failure and/or alarm	Possible Cause	Solution
Pump is not functioning	Plug not or only partially inserted in wall socket	Fully insert the plug into the wall socket
	Socket faulty or power supply out of order	Change socket or restore power supply
	On/off button not engaged	Press the on/off button
	The power cable is damaged	Replace power cable
	1 or 2 fuses out of order	Change the fuse(s)
	<b>If the power light is on and none of the above causes are found, contact your local service provider or Nausicaa Médical's After-Sales Service.</b>	
No alternating pressure	Static mode selected	Select dynamic mode
	Internal problem	<b>Contact your local service provider or Nausicaa Médical after-sales service</b>
Permanent high pressure	Internal problem	<b>Contact your local service provider or Nausicaa Médical after-sales service</b>
No or little air at pump cells	Air filter dirty Pump out of order.	Check air filter Replace pump
Low pressure in support cells	- The pressure set is too low for the patient's weight and/or position (semi-seated or seated). - Patient's weight is greater than that allowed by the pump	- Increase pressure - Adapt according to position  - Use a weight-bearing pump
Low pressure alarm	CPR open Connector air leak Cell air leak	Close CPR Check connector Check cells and connections

# Symbols (product label and/or packaging)

## 19. Symbols (product label and/or packaging)

	Indicates that the operating instructions must be read.		Indicating an important warning
	Advice or information the user and/or patient should be aware of		Indicates that the instruction manual must be read before using the product.
	Indicates that the product is a medical device compliant with EU Medical Device Regulation 2017/745		Indicates that the product complies with European Union legislation.
	Unique identification code for medical devices		Manufacturer contact details
	Production date		Serial number
	Article reference		Meets IEC 60601-1 type BF standard
<b>IP21</b>	Protection class		Indicates fuse type
	Electrical insulation class II: double insulation		Do not dispose of with household waste
	Do not use chlorine bleach (Cl)		Personal weight

# Technical features

## 20. Technical features

Model	NAUSIFLOW 2S
Reference	NA2S-COMP-PM
Class MDR	Classe I
Basic UDI	37014294NAUSIFLOW2SVX
UDI-DI	03701429400545
Electrical protection class	 Class II, Type BF, (Pas d'AP/APG)
Noise pollution	27-28 dB (maximum flow at 1 m) (ISO 7779)
Mains voltage	230 V
Sector frequency	50 Hz
Fuse rating	250 V, T1A
Protection class	IP21
Air flow	7 L/min.
Pressure range	15 to 50 mmHg
Cycle time	12 minutes
Housing material	ABS plastic
Dimensions	(L) 280 mm x (L) 205 mm x (H) 135 mm
Weight	2,6 kg

ENGLISH

## Disposal

### 21. Disposal

Follow national regulations for the disposal of any component of this pump or for its complete disposal.

## Lifespan

### 22. Lifespan

The service life indicated below depends on the frequency and intensity of use, cleaning, disinfection, transport, annual maintenance, etc.

Product	Lifespan
Pump NAUSIFLOW 2S	± five (5) years

 **The fact that Nausicaa Médical specifies a lifetime does not constitute an additional warranty.**

# Warranty

## 23. Warranty

**Article 1** : NAUSICAA Médical S.A. guarantees this device against all manufacturing and assembly defects in its components and only for devices used under the conditions provided by NAUSICAA Médical S.A. This guarantee, the conditions of which are defined below, is valid for 24 months from the date of first purchase from NAUSICAA Médical S.A. The fact that Nausicaa Médical specifies a lifetime does not constitute an additional guarantee.

**Article 2** : The warranty entitles the purchaser to free labor and free replacement of any parts found to be defective.

**Article 3** : Shipping costs to return this device and all related costs are at the reseller's expense. Goods travel at the reseller's risk and responsibility.

Under warranty: NAUSICAA Médical S.A. will bear the cost of return shipment.

Out of warranty: the reseller will bear the cost of return shipment, whether or not he accepts the repair estimate.

**Article 4** : The warranty does not apply to claims resulting from :

- Accident, misuse or negligence on the part of the purchaser.
- Transport of the device without adequate protection.
- Modifications or alterations not approved by NAUSICAA Médical S.A.
- The impact of external factors (natural disasters, fire, impact, humidity, flooding, lightning, etc.).
- Installation and/or use that does not comply with technical and safety standards, if the device is to be used in a country other than the country of purchase; and/or if the power supply is not adapted to the device's operating voltage.
- Failure to perform routine maintenance.

**Article 5** : : The reseller will not be able to invoke the benefit of the warranty:

- If the device's serial number has been removed, altered or rendered illegible.
- If the device under warranty has been modified without the approval of NAUSICAA Médical S.A.

**Article 6** : No equipment will be loaned out while defective equipment is being repaired.

**Article 7** : All warranty claims must be made through the reseller. If this is not possible, the purchaser may send the equipment directly to NAUSICAA Médical S.A. In this case, the purchaser must include the retailer's contact details and a copy of the purchase invoice in a letter enclosed with the equipment.

**Article 8** : The shipment of spare parts under warranty will only be carried out after consultation with NAUSICAA Médical S.A.'s After-Sales Service.

Please note that defective spare parts must be returned to NAUSICAA Médical S.A.'s After-Sales Service within 15 days, failing which they will be invoiced 1 month after NAUSICAA Médical S.A. has sent the parts.

**Article 9** : Defective parts replaced under or outside the warranty are guaranteed for 6 months from the date of repair or shipment of spare parts.

**Article 10** : No retailer may unilaterally modify the terms of the present warranty.

# NAUSICAA

MEDICAL

Head Office Sales Department

Phone : 04 66 51 50 80

Fax : 04 66 51 50 47

Mail : [contact@nausicaa-medical.com](mailto:contact@nausicaa-medical.com)

[www.nausicaa-medical.com](http://www.nausicaa-medical.com)

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