

GROUND PICK UP

Compatible slings

NAUSICAA Médical offers two types of slings for picking up from the ground:

- Hammock slings
- U-shape slings



When picking up the patient from the floor, it is advisable to :

- Secure the patient's neck with a rigid headrest or a headrest with fasteners.
- If the sling does not have a head support, the patient's neck should be with your hand.
- The patient's head should face away from the mast to avoid shocks during lifting.
- To make it easier for the patient's feet to pass through the lift, they should be spread out, the first along the body and the second under the bend of both knees.



DEMONSTRATION VIDEO

SCAN ME

