www.nausicaa-medical.com

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www.nausicaa-medical.com

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MATTRESS **NAUSIFLOW 512**

FOR PATIENTS FROM 30 TO 150 KG
CLASS I MEDICAL DEVICE

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Summary

User manual / Mattress : NAUSIFLOW 512

1. Introduction	32
2. Content of the box	32
3. Safety information and instructions	33
3.1 Symbols	33
3.2 Important	33
3.3 Safety instructions	33
3.4 Precautions	34
4. Target groups	34
4.1 Target patients	34
4.2 Users	35
5. Operating principle	35
6. Clinical principles	35
6.1 Objective	35
6.2 Expected clinical benefit	35
6.3 Indications	36
6.4 Contraindications	36
6.5 Other measures required	36
6.6 Effects	36
6.7 Manual control	36
7. Presentation	37
7.1 Composition	37
7.2 Description	38
8. Installation	38
8.1 Precautions & restrictions on use	39
9. Uninstalling	39
10. Maintenance	40
10.1 Daily & weekly maintenance	40
10.1.1 Daily procedure	40
10.1.1.1 Cover maintenance	40
10.1.1.2 Cover disinfection	40
10.1.2 Weekly procedure	41
10.1.2.1 Care of the cover and mattress	41
10.1.2.2 Mattress disinfection	41

11. Pre-storage maintenance	42
11.1 Generality	42
11.2 Check	42
11.2.1 Air mattress & cover	42
11.3 Cleaning procedure for the cover and mattress	42
11.4 Mattress disinfection	42
12. Storage	44
13. Climatic requirements	44
14. Compatibles pumps	44
15. Troubleshooting	45
16. Symbols (product label and/or packaging)	46
17. Technical data	47
17.1 Mattress	47
17.1.1	47
17.2 Foam mattress underlay	47
17.2.1	48
17.3 Cover	48
17.3.1	48
18. Elimination	48
19. Lifespan	49
20. Warranty	50

Introduction

1. Introduction

Dear customer,

Nausicaa Medical would like to thank you for choosing this mattress.

To avoid any damage or danger caused by incorrect use, please read this manual in its entirety before using the mattress for the first time and before using it again. This manual contains important information and notes necessary for the correct use of the mattress. If you have any questions about safety warnings, precautions or use, contact your local service provider or Nausicaa Médical's After Sales Service. To avoid injury or damage due to misuse, do not use the mattress until all your questions have been answered.

These instructions must always accompany and/or be supplied to any person requiring the use of this mattress or to the person responsible for using it.

The life of this mattress depends essentially on correct installation and use.

This product complies with the requirements of the following standard(s):

EN 60601-1 : 2006/A1 : 2013

EN 60601-1-2 : 2015

EN 60601-1-11 : 2015

Regulation (EU) 2017/745 of 05 April 2017

 The term «system» used in this manual refers to all the components, i.e. « this mattress combined with a compressor ».

Content of the box

2. Contenu de l'emballage

Please check that the content listed below corresponds to the content you received.

1 13 cm air mattress + foam underlay

1 user manual 

Safety information and instructions

3. Safety information and instructions

3.1 Symbols

To make this user manual easier to read and understand, here is a brief explanation of the main symbols used.

Some of these symbols and their explanations require special attention.

 **USER'S MANUAL:**
You must read this manual before using the machine.

 **INFO:**
Provides advice or information of which the user and/or patient should be aware.

 **ATTENTION/DANGER :**
Indicates an appropriate use, operation or maintenance procedure to avoid altering and/or rendering inoperative this mattress or one of its components or damaging other property.
Draws attention to a potential hazard requiring appropriate procedures or practices to prevent injury.

For more information on the symbols used, see page 46.

3.2 Important

Read and follow the safety instructions carefully. They consist of a text or a combination of a symbol and a text. The symbol used does not replace the text!

The safety instructions apply to anyone using and/or intervening in any way whatsoever on this mattress or one of its components.

The person who makes this mattress available to the user and/or a patient must inform them of the potential risks and train them in its correct use so that no additional risk results from inappropriate behaviour.

 Take particular care when children, disabled people or pets are in the vicinity of the mattress. Use this mattress only for its intended purpose. Any other use is prohibited.
DO NOT use with compressors other than those compatible with and marketed by Nausicaa Medical.

3.3 Safety instructions

 Repairs and/or work on this mattress may only be carried out by qualified personnel or by personnel who have received the necessary training. No modifications are permitted and only original parts may be used.

Do not use this mattress near flammable liquids or gases. Keep it away from all sources of heat (heaters, electric blankets,) and naked flames (candles, etc.).

To prevent damage to the compressor power cable, insert it into the cable duct provided.

Safety information and instructions

To maintain the mattress's effectiveness, do not add any extra thickness between the patient and the mattress. Simply covering the mattress with a sheet without tucking it in is more than enough.

Check that the bed rails are appropriate for the correct use of this mattress by ensuring that the space between the top of the uncompressed mattress and the top edge of the top rail is at least 22 cm and that the space between the rails is as small as possible to avoid any entrapment zones.

Failure to do so could result in serious injury or even suffocation.

Do not use any components other than those recommended by Nausicaa Médical. Use of a non-compliant component will automatically invalidate the legal warranty.

i It is recommended that you notify the service provider, the manufacturer and the competent authority of any serious incident* that occurs during the use of this mattress or if you consider or have reason to believe that it is or has been tampered with.

*Serious incident» means any incident which has directly or indirectly led to, or could lead to:

- (a) The death of a patient, user or other person;
- (b) the serious deterioration, temporary or permanent, of the state of health of a patient, user or other person;

3.4 Precautions

- Protect this mattress from sharp objects and surfaces, including the claws and teeth of animals.
- Before use or re-use, this mattress must be clean and dry.
- Follow the recommendations for care, cleaning and disinfection (pages 40 to 44).

i Store this mattress as indicated in this user manual (page 44).

i Observe the minimum and maximum weights that this mattress can support. Failure to do so may cause serious harm to the patient, lead to rapid deterioration, cause damage and render the mattress inoperative.

Target groups

4. Target groups

4.1 Target patients

The target patients are mainly adults with a height ≥ 150 cm and a weight between 30 kg and 150 kg who, because of illness, injury, disability or age, are bedridden for ≥ 15 hours a day, have one or more pressure sores or are considered to be at risk because of a temporary or permanent deterioration in their state of health.

Target groups

4.2 Users

Users are generally healthcare professionals, who may be assisted by carers such as care assistants or family members.

That's why this mattress is designed to be simple and practical to use and maintain, whoever the user may be.

For use by a member of the family, it is essential that they have received full and appropriate instructions from experienced staff or the installer, and that they are able to carry out the supervisory role.

Operating principle

5. Operating principle

The mattress cells are inflated alternately according to a specific pattern governed by the associated compressor.

This pattern creates a cyclical and alternating pressure movement on the different areas of the body.

Clinical Principles

6. Clinical principles

This mattress is one of the components of a system intended to be used for medical purposes to prevent the appearance of a pressure sore and/or to provide assistance in caring for it.

6.1 Objective

Reduce pressure on the areas of the body most prone to developing a pressure sore.

6.2 Expected clinical benefit

Maintains tissue oxygenation in anatomical areas in contact with the mattress by reducing the pressure applied to cutaneous and subcutaneous tissues.

6.3 Indications

This mattress is mainly recommended for the prevention and/or treatment of stage 1 to 2¹ pressure sores.

Its aim is to reduce the incidence of pressure sores while optimising patient comfort. It is particularly suitable for patients who are bedridden for ≥ 15 hours a day and get up during the day, and whose risk is assessed as 'medium to high*'. In addition to care, it is indicated for stage 1 to 2¹ pressure sores.

These recommendations are based on clinical judgement and an individual protocol including position change, nutritional support, skin care, wound care, etc. All other aspects of care must be undertaken by a doctor and an experienced carer.

¹ According to medical opinion - Classification of stages according to the Braden scale.

- i** Regularly monitor the compatibility between the system and the patient's state of health, particularly if their physical condition has changed (amputation, weight gain/loss, mental changes, etc.).
Always take a holistic view of the patient's state of health.

Always take into account any contraindications. If necessary, consider using a mattress type 834 or 1034, respectively 20 and 25 cm all-air.

6.4 Contraindications

Do not use the device if the patient has one or more of the following contraindications:

- Less than 30 kg or more than 150 kg,
 - non-stabilised fracture of the spinal column,
 - cervical or skeletal traction,
 - multiple acute trauma,
 - non-stabilised post-traumatic fracture
- If this is the case, use is subject to the advice of a qualified doctor.

- i** It is also recommended not to use this system in a hyperbaric chamber or on a stretcher.

6.5 Other measures required

The system alone is not enough; other preventive measures are also essential, such as:


- Changing position at least every 2 to 3 hours;
- Maintaining good skin hygiene and avoiding or reducing maceration as much as possible;
- Changing protection regularly, especially in the case of severe incontinence;
- Observing the state of the skin every day or having someone else do so;
- Ensure that the patient drinks regularly and in sufficient quantities.

6.6 Effects

There are no scientifically-documented adverse effects, although there may be a slight risk that the alternating action may cause and/or promote spasticity.
In such cases, the advice of a specialist doctor should be sought.

6.7 Manual control

In alternating mode, to check that the pressure is correctly adjusted, you should locate the least inflated cell under the patient's buttocks and try to insert a hand flat palm up (see illustration below).

6.7 Manual control

You need to be able to feel slight pressure while still being able to move back and forth. The acceptable range between the two is between 25 and 40 mm. This manual checking procedure is issued by the Agency for Health Care Policy and Research (AHCPR).

-  If you are completely unable to move back and forth, this means that the pressure is too low (heeling) > **increase the pressure.**
-  If you can perform this movement without any resistance, or even with extreme ease, this means that the pressure is too high < **reduce the pressure.**

Presentation

7. Presentation

This mattress can be used in short- and long-term care settings, both in-patient and out-patient, in retirement homes, nursing homes, medical-social establishments, etc., as well as at home.

7.1 Composition

This mattress consists of 17 independent, disconnectable nylon/PU cells with a total air height of 13 cm, resting on a 5 cm thick foam mattress pad.

For optimal comfort, the three cells located in the head area are static, and to facilitate pressure relief in the heel area, the last four cells are equipped with an automatic quick release closure connector.

Its air supply tubes are pinch-proof, and its compressor connection is equipped with a transport cap.

In transport mode (transport cap on the connector), its static air autonomy is ± 6 hours*. It is equipped with a side cable range and a **CPR** valve that allows rapid deflation (< 20 seconds) in case of extreme emergency (cardiac arrest).

This is located in the head area on the left side* (*viewed from the foot of the bed).

Its cover is made of durable Poly MeSH with a non-slip underside. It has two straps for securing the mattress to the bed base.

It is covered with a fully zipped bi-elastic polyurethane cover. It is liquid-proof, breathable, treated with silver ions (Ag+), classified as M1 fire-resistant and machine washable at a maximum temperature of 95°C.

-  In the event of a power failure, the first recommendation is to switch to transport mode (see above). However, the mattress can remain connected to the compressor during the power failure as long as the compressor was not in the process of changing cycles.
WARNING! Autonomy calculated based on a patient weighing between 60 and 80 kg and measuring between 160 and 180 cm, provided that the bed base is completely flat.

Presentation

7.2 Description

Cellules NPU

Connector +
transport cap

Self-sealing quick
connector

CPR valve

Fastening strap

Cable conduit

Zipped cover

Installation

8. Installation

Before installation, always check that the mattress is not damaged.

DO NOT install or use a mattress that is even slightly damaged!

- Remove the existing mattress from the bed.
- Place the air mattress directly on the bed base.
- Unroll the mattress on the bed base, ensuring that the air supply connector is at the foot of the bed.
- Using the two straps, secure the mattress to the bed's backrest.
- Check that the mattress is properly inserted between the handles located on the outside of the bed base frame.
- Check that all moving parts of the bed are working properly.
- Check that the CPR valve on the mattress is tightly closed.
- If available, place the compressor's electrical cable in the cable duct provided so that it does not pose a hazard and cannot be caught, crushed, trapped or pulled out.
- Insert the mattress air supply hose connector into the compressor connector.
- Allow the mattress to inflate completely¹ before placing the patient on it.
- Cover the mattress cover with a sheet that is as thin as possible.

⚠ Ensure that both connectors (mattress and compressor) are properly connected and secured.

⚠ Tucking the sheet too tightly reduces the effectiveness of the mattress.

Installation

- i** ¹ Please note that the inflation time depends on the dimensions of the mattress (L x W x H).
Example: it takes approximately 25 minutes for an 88x13*x200 cm mattress to be fully operational when using the Nausiflow 100 AUTO QUATTRO compressor.
*Air height

Depending on the mattress used, the Nausicaa rapid inflation pump will allow you to have it ready for use in ± 2 minutes.

ENGLISH

8.1 Precautions & restrictions on use

Take all necessary precautions

- When the mattress is used at home without the assistance of healthcare professionals or similar personnel (caregivers, etc.).
- When used during the first few days after surgery for bedsores (skin graft or flap). **In this specific case, it is preferable to use the static mode available on all our compressors.**

Uninstalling

9. Uninstalling

- Switch off the compressor.
- Remove the power cable plug from the wall socket.
- Remove the power cable from the cable duct.
- Disconnect the mattress from the compressor.
- Open the **CPR** valve and allow the mattress to deflate².
- Once the mattress is deflated, remove the fastening straps, place the compressor in the centre of the mattress and carefully roll up the mattress.
- Place the mattress and compressor back into the carry bag.

² The pump offered by Nausicaa Médical also allows you to deflate the mattress in ± 2 minutes.

10. Maintenance

It is MANDATORY to clean and disinfect the mattress in accordance with applicable national and/or local guidelines for cleaning and disinfecting medical devices. In this section, we describe the procedures to be followed. It is important to follow them.

-  The cover must be cleaned daily using a clean cloth soaked in a cleaning/disinfectant solution. It is not necessary to disconnect the mattress from the compressor during this operation.
-  To maintain the maximum hygiene of the mattress cover, it must be machine washed at least once a week.
-  To maintain maximum hygiene, the mattress must be cleaned thoroughly (inside and outside) at least once a week.

As mattresses can harbour germs, we recommend that you wear protective gloves **BEFORE** handling, cleaning or disinfecting them.

10.1 Daily & weekly maintenance

-  Before handling and/or working on the mattress, put on protective gloves and disinfect the surface.
-  If you notice a hole, wear or any other damage to the cover or any of the mattress components, replace it immediately.

10.1.1 Daily procedure

10.1.1.1 Cover maintenance

If possible, remove the patient from the mattress and proceed as follows:

Check the integrity of the cover.

- Clean any traces of faeces and blood first.
- Soak a clean, single-use wipe or cloth in cleaning solution and scrub the entire surface, paying particular attention to the area where the patient is lying. Repeat if necessary.
- Remove any excess cleaning solution with a clean, slightly damp single-use cloth.
- Wipe with a clean, dry single-use cloth.

10.1.1.2 Cover disinfection

- Spray a disinfectant solution over the entire surface of the cover and leave it to act according to the manufacturer's recommended instructions for use.
- After the recommended waiting time, wipe off any excess with a clean, dry, single-use cloth if necessary.
- Ensure that there is no residual moisture or condensation and that the cover is completely dry before repositioning the patient on the mattress.

Maintenance

If you are unable to remove the patient from the mattress, proceed as follows:

- Check the integrity of the cover as much as possible.
- Turn the patient onto their side, making sure to secure them with bed rails and, if necessary, prevent them from rolling back onto their back using a device (positioning cushion or other).
- Clean any traces of faeces and blood first.
- Soak a clean, single-use wipe or cloth in a cleaning/disinfectant solution. Scrub the accessible surface, paying particular attention to the area where the patient is lying, and leave to act according to the manufacturer's recommendations for use.
- Ensure that the patient, and in particular their skin, does not come into contact with the cleaning/disinfectant solution during the entire operation.
- Wipe off any excess with a clean, dry, single-use cloth.
- Repeat the entire process above for the other side.
- Ensure that there is no residual moisture or condensation and that the cover is completely dry before repositioning the patient on the mattress.

10.1.2 Weekly procedure

10.1.2.1 Care of the cover and mattress

Remove the patient from the mattress and proceed as follows:

- Turn off the compressor (power indicator light off), disconnect the mattress from the compressor and place the transport cap on the connector throughout the operation.

Cover :

- Remove the cover from the mattress completely and check that it is intact.
- Machine wash using a cleaning/disinfectant solution suitable for this type of care, ensuring that the maximum washing temperature does not exceed **95°C**.

Mattress :

- Check the integrity of the mattress.
- Soak a clean, single-use wipe or cloth in a cleaning solution and rub all its internal components (grey base, cells, tubes and CPR valve).
- Repeat this operation on all its external components (grey base, air connector and CPR strap).
- Remove any excess cleaning solution with a clean, slightly damp disposable cloth.
- Wipe with a clean, dry disposable cloth.

10.1.2.2 Mattress disinfection

- Check the integrity of the mattress.
- Spray a disinfectant solution on all its internal components (the grey base, cells, tubes and CPR valve) and then use a sponge or clean disposable cloth to scrub all these components. Leave the solution to act according to the manufacturer's instructions for use.
- Repeat the same process as above on all external components (base, air connector and CPR strap).
- After the recommended exposure time, if necessary, wipe off any excess with a clean, dry, single-use cloth.
- Allow to dry, then ensure that there are no traces of residual moisture or condensation and that the mattress is completely dry before covering it with its cover.

Pre-storage maintenance

11. Pre-storage maintenance

11.1 Generality

It is recommended that you check and replace any defective parts and perform a complete maintenance check before storing the mattress. This work must be carried out by qualified personnel.

11.2 Checks

11.2.1 Air mattress & cover

Check the following items for visible damage and/or signs of significant wear:

- The cells, their connectors and fastenings (press studs and translucent rings)
- The tubes
- The non-return valves
- The CPR valve and strap.
- The transport mode cap and connector
- The base (grey cover) of the mattress
- The waterproof compartment for the foam under-mattress and its zip
- The condition of the foam
- The cable sleeve and its press studs
- The fastening/rolling straps
- The removable cover and its zip

If the mattress is in good condition, move it to the cleaning area and place it on a stable, firm work surface.

11.3 Cleaning procedure for the cover and mattress

- Remove the foam mattress pad from its compartment.
- Completely remove the mattress cover.
- First remove any traces of faeces and blood, then machine wash using a cleaning/disinfectant product suitable for this type of care, ensuring that the maximum washing temperature does not exceed 95°C.
- Lay the mattress flat or hang it up, then spray a cleaning solution on all its internal components (grey base, cells, tubes, CPR valve and inside the compartment housing the foam mattress pad) and scrub all these components with a sponge or clean disposable cloth.
- Repeat the same process as above for the outside of the mattress (base, air connector and CPR strap).
- Remove any excess cleaning solution with a slightly damp clean, single-use cloth.
- Wipe with a clean, dry, single-use cloth.

11.4 Mattress disinfection

- Spray a disinfectant solution on all internal components (grey base, cells, tubes, CPR valve and inside the compartment housing the foam under-mattress). Leave the solution to act for a few moments (see manufacturer's recommendations for use) then, using a sponge or clean disposable cloth, scrub all these components.
- Repeat the same process as above on all external components (base, air connector and CPR strap).

Pre-storage maintenance

- After the remanence period, if necessary, wipe off any excess with a clean, dry, single-use cloth.
- Store the mattress in a suitable place to eliminate all traces of residual moisture or condensation.
- Make sure that the mattress is completely dry before covering it with its cover and putting the foam underlay back into its compartment.
- Once the mattress has been reconditioned, check that everything is correct, then connect it to a compressor and inflate it to maximum (Care Mode) to make sure there are no leaks.
- After checking that everything is working correctly, disconnect it from the compressor, open the CPR and allow the mattress to deflate before rolling it up, placing it in a clean plastic bag and inserting it in its carry bag.
- Affix an identification label to the transport bag to ensure that the mattress is cleaned, disinfected and in good working order.
- Store the transport bag as recommended.

 You can use a cleaning/disinfecting solution, but make sure that the disinfecting agents have a sufficiently wide range of action.

 **PRECAUTION:** when it comes to disinfecting the cover, the polyurethane (PU) coating on our covers is generally highly resistant to alcohol-based disinfectants. Phenol or its derivatives can also be used in dilute concentrations (follow the recommendations of the manufacturer of the product used). However, we always advise you to test the product before use.
High concentrations of active chlorine can have a negative effect on PU-coated materials. Its use is therefore not recommended.
If you decide to use it, use very dilute concentrations and rinse thoroughly with water. The same applies to peroxide-based disinfectants and, in general, to all disinfectants based on an oxidising reaction mechanism.

 **DRYING:** Proper drying is a crucial factor in the care and maintenance of coated items. This can be done by hanging for the base and its components and by machine at moderate temperatures for the cover.

 If machine washing is necessary, the maximum temperature should not exceed 95°C.

 The use of fabric softeners is not permitted.

 As the cells in the heel area are fitted with a quick connector, ensure that they are securely fastened to prevent liquid from entering the tubes and cells.
The same applies to the other cells and the CPR valve.

 During daily or weekly maintenance, do not press on the inside of the flap protecting the zip of the compartment housing the foam mattress underlay, or on the zip itself, to prevent liquid from entering this compartment.
If this happens, the foam underlay must be changed and the new one inserted only after checking that there is no residual moisture or condensation inside the compartment and that it is completely dry.

 Sterilisation can have a negative effect on polyurethane (PU) coverings. It is therefore not recommended.
If you opt for this treatment, bear in mind that the product's lifespan will be shorter. This is unavoidable for all thermoplastic polyurethanes, due to their chemical composition.

Storage

12. Storage

We recommend :

- Store the mattress in a carry bag or in the original packaging.
- Handle the device with care.
- Do not place heavy objects on it.
- Protect it from direct sunlight.
- Store in the correct climatic conditions.

ENGLISH

Climatic requirements

13. Climatic requirements

Conditions	Temperatures (°C)	Humidity level (%)	Atmospheric pressure (hPa)
Use	From 10°C to 40°C	From 30% to 75% (non-condensing)	From 700 hPa to 1.060 hPa
Storage/transport	From 0°C to 40°C		

Compatibles pumps

14. Compatibles pumps

Mattress di- mensions in cm (LxW)	Model	Maximum load weight
198 x 70	Nausiflow 2S Nausiflow 100 AUTO QUATTRO R'GO SOINS Nausiflow 100 MAXI 2	150 kg
198 x 88		
198 x 100	Nausiflow 100 AUTO QUATTRO R'GO SOINS Nausiflow 100 MAXI 2	150 kg
198 x 120	R'GO SOINS Nausiflow 100 MAXI 2	150 kg

Troubleshooting

15. Troubleshooting

The following table lists some problems, possible causes, and the solution to be provided. If the solutions in the panel do not solve the problem, contact your distributor or Nausicaa Medical's after-sales service.

Failure and/or Alarm	Possible cases	Solutions
Low pressure in the cells of the mattress It is not inflated or is insufficiently inflated (the mattress no longer offers support to the patient)	The CPR is open or not hermetically sealed Or faulty.	Seal the CPR Change the CPR
	One or more cells are not correctly connected to the hoses	Connect all the cells to the respective tubes
	The 3-way connector on the mattress is not fully inserted into the hoses	Push the mattress connector until it touches the tubes. When doing so, make sure that the anti-torsion spring is not positioned askew
	The spacers on the hoses are not correctly connected.	Connect all the pipe spacers securely
	The mattress connector is not correctly connected to the compressor connector	Refasten the mattress connector to the compressor connector. A click will be heard when it is correctly locked
	One or more hoses are bent, crushed, twisted or cut	Replace the defective hoses
	There is a leak in one or more cells.	Replace the cell(s)
Muffled noise with a vibrating effect when inflating, seeming to come from the compressor	Problems with one or more of the non-return valves (yellow and red) located on the air hoses inside the mattress.	Change the valve(s)

Symbols (product label and/or packaging)

16. Symbols (product label and/or packaging)

	Instructions for use available		Indicates an important warning
	Advice or information of which the user and/or patient should be aware		Indicates that the instructions for use must be read before using the product
	Indicates that the product is a medical device compliant with the EU Medical Device Regulation 2017/745		Indicates that the product complies with European Union legislation
	Unique identification code for medical devices		Manufacturer's contact details
	Date of manufacture		Serial number
	Item reference		Do not dispose of with household waste
	Orientation marker		Cover washing - Maximum washing temperature (95 °C)
	Hand washing recommended		Mattress and cover - Do not dry clean
	Do not use chlorine (Cl) agent		Mattress - Do not tumble dry
	Person size		Mattress and cover - Do not iron
			Person weight

Technical data

17. Technical data

17.1 Mattress

NAUSIFLOW 512	
MDR class	Class I
IUD-DI basic	37014294NAUSIFLOWAUTOQGN
Cell material	Nylon Polyurethane (N/PU) - 0.2mm
Number of cells	17
Number of cells with quick connector	4 (Cells in the heel area)
Air height	13 cm
CPR deflation time	< 20 seconds
Anti-equine cushion	Not available
Base material	Poly MeSH
Base thickness	0.97mm

ENGLISH

17.1.1

NAUSIFLOW 512			
Dimensions in cm	Reference	IUD-DI	Weight including foam
198 x 70 x 18	NA512NPU-MAT70	03701429409128	6,8 kg
198 x 88 x 18	NA512NPU-MAT88	03701429409111	7,75 kg
198 x 100 x 18	NA512NPU-MAT100	03701429410971	7,95 kg
198 x 120 x 18	NA512NPU-MAT120	03701429409135	9,60 kg

17.2 Foam mattress underlay

NAUSIFLOW 512 foam underlay	
MDR class	Class I
IUD-DI basic	37014294NAUSIFLOWAUTOQGN
Type of foam	HD Polyether
Thickness in cm	5 cm
Density	20 kg/m ³

Technical data

17.2.1

NAUSIFLOW 512 foam underlay		
Dimensions in cm	Reference	IUD-DI
198 x 70 x 5	MOU-0198-01	03701429410957
198 x 85 x 5	MOU-0198-02	03701429409074
198 x 100 x 5	MOU-0198-04	03701429411053
198 x 120 x 5	MOU-0198-03	03701429410964

17.3 Cover

Mattress cover NAUSIFLOW 512	
MDR class	Classe I
Material	Polyurethane-Polyester
Thickness	0,5 mm

17.3.1

Mattress cover NAUSIFLOW 512		
Dimensions in cm	Reference	IUD-DI
200 x 70	NAUS-M51270NPU-HOU	03701429408756
200 x 88	NAUS-M51288-HOU	03701429408718
200 x 100	NAUS-M512100-HOU	03701429411060
200 x 120	NAUS-M512120NPU-HO	03701429410063

Elimination

18. Elimination

Comply with national regulations for the disposal of any of the components of this device or for the disposal of the entire device.

Lifespan

19. Lifespan

The service life indicated below depends on the frequency and intensity of use, the frequency and quality of maintenance and disinfection, as well as the products used for these two processes, transport, handling and storage.

Product	Lifespan
Nausiflow 512 mattress	± five (5) years
Foam	one (1) year
Cover	± two (2) years

 The fact that Nausicaa Médical specifies a lifetime does not constitute an additional guarantee.

Warranty

Article 1 : NAUSICAA Medical S.A. warrants this device against any manufacturing defects and assembly of its components, and only for devices used under the conditions provided by NAUSICAA Medical S.A. This warranty, the conditions of which are defined below, is valid for 24 months from the date of the first purchase from NAUSICAA Medical S.A.

Article 2 : The warranty entitles you to free labour and free replacement of parts found to be defective.

Article 3 : The outbound port of the device, as well as all related costs, are the responsibility of the reseller. The goods always travel at the risk and under the responsibility of the dealer.

Sous garantie : the costs of return after intervention will be borne by the company NAUSICAA Medical S.A.

Hors garantie : the return costs are at the expense of the retailer whether or not he accepts the repair estimate.

Article 4 : The warranty does not apply if the claims are resulting from:

- Accident, misuse of the appliance, or negligence of the purchaser.
- Carrying the device without adequate protection.
- A modification or transformation not validated by the company NAUSICAA Medical S.A.
- The impact of external agents (natural disaster, fire, shocks, humidity, floods, lightning, etc.).
- Installation and/or use in a manner that does not comply with technical and safety standards if the device is to operate in a country other than the country of purchase and/or if the power supply is not suitable for the operating voltage of the device.
- A lack of routine maintenance.

Article 5 : The reseller will not be able to invoke the benefit of the warranty:

- If the device's serial number has been removed, altered, or rendered illegible.
- If the device under warranty has been modified without the approval of NAUSICAA Medical S.A.

Article 6 : During the repair of defective equipment, no equipment loan will be made.

Article 7 : Any warranty claim shall be exercised through the reseller. If this is not possible, the buyer may send his equipment directly to the company NAUSICAA Medical S.A. In this case, the buyer must indicate in a letter attached to the equipment, the contact details of the dealer and a copy of the purchase invoice.

Article 8 : The shipment of spare parts under warranty will be made or will only be made after consultation with the After-Sales Service of NAUSICAA Medical S.A.

Please note that defective spare parts must be returned to the After-Sales Service of NAUSICAA Medical S.A. within 15 days, otherwise they will be invoiced 1 month after the parts have been sent by NAUSICAA Medical S.A.

Article 9 : Defective parts changed under or out of warranty are guaranteed for 6 months from the date of repair or shipment of spare parts.

Article 10 : No dealer may unilaterally change the terms of this warranty.

NAUSICAA

MEDICAL

Sales Department Headquarters

Phone : 04 66 51 50 80

Fax : 04 66 51 50 47

Mail : contact@nausicaa-medical.com

www.nausicaa-medical.com

Your contact :