



www.nausicaa-medical.com

This medical device is intended for patient holding and positioning. It is not a restraint system. Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

SLIM FIT

FOR ALL CHAIRS




TO PREVENT FORWARD
SLIDING AND MAINTAIN THE
CHEST IN UPRIGHT POSITION
FOR ALL CHAIRS EXCEPT
SHELLS



Maximum useful area: + 230 cm
Belt height and width
pelvic part: + 11 cm
Width strap: + 48 mm
Buckle: 50 mm

Its shoulder heights are separately adjustable to help position and compensate for a bad position.

- Available in 4 sizes S, M, L, XL
- Maximum useful area: 230 cm (the backrest and seat of the chair must have a separation space).
- Waist (in cm)



The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative. This medical device must be implemented or trained by a health professional.



CSYP2HP-PF

v 05.02/23

Important Recommendations :

- In order to get a maximum efficiency of these sort of product, it is necessary :
- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient
- These products must not be in direct contact with a wounded skin.



SLIM FIT

HOSPITAL RANGE

Combination of polyurethan foam
and on both faces
Jersey/PU coating
100% polyester bias
100% polyester strap
Nylon buckle
Thermo compressed foam

- Position the pelvic part on the seat of the chair and let the buckle go out towards the back of the chair.
- Let the patient sit on the chair.
- Position the belt around the patient's waist.
- Close the belt using the velcro strap.
- Place the chest belt being careful that the shoulders pads are positioned on the patient's shoulders.
- Pass the strap of the back part of the belt behind the backrest of the chair.
- Before closing it and without adjusting it, check that the back strap has a female buckle looking downwards and 2 female buckles looking upwards.
- Close the 4 rear buckles, adjust