NAUSICAA

MEDICAL

www.nausicaa-medical.com

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www.nausicaa-medical.com

ENGLISH

STAND UP LIFT: **EASYLEV COMPACT** OR WITH MECHANICAL OR ELECTRIC ECP

EN

USER MANUAL

STAND UP LIFT: **EASYLEV FIXED OR WITH MECHANICAL OR
ELECTRICAL OPENING LEGS** ————— p 31-55

Summary

Assembly Guide / Stand Up Lift: EASYLEV

• Presentation, contraindications, contents	37
• Setup instructions	38-40
• General safety instructions	41
• Safety instructions before use	42-43
• Use of the device	44-45
• Use of the motorization	46-47
• Use of the sling	48-50
• Storage	51-52
• Reuse of the stand up lift	53
• Labeling	54-55
• Cleaning and maintenance	56
• Preventive maintenance and safety check	57-60
• Technical and dimensional characteristics	61
• Spare parts	62
• Motorization technical characteristics	63-65
• Troubleshooting guide	66
• Warranty	67

The Easylev stand up lift is a class I medical device in accordance with Regulation (EU) 2017/745.

ALL OUR DEVICES COMPLY WITH THE NF EN ISO 10535:2021 STANDARD.

Presentation

The EASYLEV stand up lift is a class I medical device in accordance with Regulation (EU) 2017/745.

It can be used in hospitals, clinics, nursing homes, care homes, at home... and must be operated by a caregiver.

It is designed for use with an appropriate sling.

This device has been designed for patients weighing up to 160 kg.

The climatic conditions for using the EASYLEV are as follows: ambient temperature from 0°C to 40°C, humidity from 20% to 80%, and air pressure from 700 to 1060 hPa.

The device can be used in the bathroom (without contact with water) or in toilets.

The stand up lift allows the patient to be transferred from a bed, a recliner, a toilet, or a wheelchair.

Contraindications

Diseases such as osteogenesis imperfecta, osteoporosis, spinal disease, mental illness, or epilepsy may be contraindications.

Contents

The Easylev stand up lift has already been inspected to ensure it is free of defects and that nothing is missing. Nevertheless, check the product immediately upon receipt for any damage that may have occurred during transport. Use the delivery note to check that all items are present and that the delivery is therefore complete.

Product concerned	Box contents	Quantity
EASYLEV	Stand up lift	1
	Control unit	1
	Battery	1
	Remote control	1
	Knee support	1
	QR code (user manual)	1

Setup instructions

BEFORE USING YOUR STANDER, IT IS NECESSARY TO CHECK THAT:

- The wheels turn and roll normally.
- The operation of the rear wheel brakes is correct.
- There is no wear or deformation on the hooks.
- All the different components are properly connected.
- The battery is sufficiently charged.
- The stander's legs open and close properly (only for the "ECP" version).

WARNING:

- Unfolding must be done with the brakes locked for the "Compact" and "ECP" models and with the legs tightened only for the "ECP" model.
- Do not trap the remote control cable between two metal parts.



The device is delivered folded in half.

Step 1: Loosen the upper knobs to release the mast **without removing them.**

Step 2: Completely unscrew the lower knobs and **remove the mast axle.**

Step 3: Remove the strap from the knee support wrapping around the base of the device, then raise the upper part of the device.

Setup instructions


Step 4: Put the axle back in place as well as the washers and knobs in the lower hole.

Step 5: Tighten the four knobs firmly to lock the assembly.

Step 6: Slightly loosen the knob then remove the knee support from its slot.

Step 7: Insert the knee support into its housing and tighten its knob to secure it, then lower the mast.

ENGLISH

The device is ready to be used.

Setup instructions

To fold the device, apply this procedure in reverse order.

User instructions:

- Your stander is designed to transfer a patient. it should not be used for other purposes.
- Check that the patient's weight does not exceed the maximum weight the stander can support.
- Manoeuvre the stander using the handlebar, never by pushing the patient.
- The stander must be handled with care when transferring a patient and at a speed appropriate to the situation.
- Move the stander on flat and smooth surfaces. it is not recommended to use it on a slope of more than 5°: if you must use a ramp, it is advisable to have a second person assist you. never use a stander in a shower.
- Never recharge the battery of a stander near a bath or shower.
- It can be used in damp rooms such as bathrooms and toilets. it is not designed to be used in a shower.
- Standeres are useful and effective when a patient has a certain degree of dependency and needs a caregiver to perform basic care.
- They contribute to:
 - stimulating the patient and their mobility;
 - stimulating the cardiac system;
 - combating osteoporosis and disorders associated with immobility such as falls during transfers;
 - stimulating brain activity, thus promoting the maintenance of urinary and fecal continence;
 - improving bowel activity and bladder function.

General safety instructions

1. Use the stander only for its intended purpose, in compliance with current legislation for medical devices, the legal requirements concerning them, regulations for occupational safety, and accident prevention.
2. Note that the stander is a medical device, so its user is required to comply with the regulations on the use of medical devices.
3. The requirements for the electrical installation of the room or area where the stander is used must comply with current regulations.
4. Use the stander only after having been trained in its handling and with full knowledge of the facts.
5. Before commissioning, read the complete user manual to avoid damage due to improper handling or exposure to risks. The manual contains important information and remarks necessary for the use of the stander.
6. Use the stander only in accordance with this manual. It can be consulted at any time via the QR code on the mast. Keep the manual carefully for future reference in case of questions. Attach this manual to the stander in case of a change of owner.
7. Before each use, check that the stander and its accessories are in working order and in perfect condition.
8. Before using the stander with other medical or non-medical devices, check that the combination of these products is authorized and that they can be used together safely.
9. Assembly, commissioning, maintenance, and repair of the stander must be entrusted only to competent specialists.
10. It is the responsibility of the user or operator to ensure (by appropriate measures and instructions) that electrical stress on the charging cord (by bending, pulling, shearing, crushing) is excluded during charging or cleaning of the room.
11. Observe the operating time and the maximum permissible load. These values must not be exceeded, otherwise safe operation can no longer be guaranteed.
12. Do not expose the stander to direct sunlight, heat, or humidity.
13. Ensure that no moisture enters the electrical system.
14. Avoid electrical stress on the electrical cords used. Pulling, bending, or crushing the electrical cords may damage them.
15. Recharge the batteries in a well-ventilated area.
16. Electromagnetic or other disturbances between the stander and other devices cannot be ruled out. If there is a risk of such disturbances, the source of the disturbances must be moved away or the stander should not be used.
17. Disturbances caused by the use of portable communication devices cannot be completely ruled out. Therefore, a safety distance of at least 3.3 m must be maintained to ensure the safe operation of the stander.
18. Do not allow children to stand unsupervised near the stand assist device.
19. The stand assist device must no longer be used if abnormal noises, damage, or any other malfunction occurs. In this case, do not connect the stand assist device to the charger, but inform NAUSICAA MEDICAL.
20. If it is damaged or defective, the stand assist device must no longer be used and must not be plugged into the mains. Inform the reseller to have them remedy the defect or breakdown.
21. Any serious incident related to the device should be reported to NAUSICAA MEDICAL and to the competent authority of the Member State in which the user and/or patient is established.

Safety instructions before use


Follow the safety instructions below before each patient transfer:

- The caregiver must have the necessary knowledge to select and use an appropriate lifting sling.
- Check the compatibility between the lifting arm hooks and the strap before using the stander.
- Check the maximum load of the device and the sling.
- Check the size and model of the sling according to the patient's body shape.
- Check the condition of the sling before each use. The sling must not have any tears in the fabric or damaged seams.
Check if the correct attachments have been fastened. All sling loops have 3 different levels: long (green) - medium (yellow) - short (red). Each pair of loops must only have the following attachment combination: long (green), medium (yellow), short (red).
- Check that all loops are attached to the hook of the lifting arm.
- Lock the wheels of the wheelchair or medical bed so you can safely lift and lower the patient. Leave the wheels of the patient's stand assist device in the unlocked position.
- Observe the patient's behavior during the transfer. Sudden movements by the patient or obstacles may cause hazards.
- Raise the patient to the required height.
- Keep the lifting straps away from intense heat or open flames; they are not fireproof.
- Before each use, it is important to check the condition of the entire transfer device (stand assist device + sling): stitching, condition of the straps and fabric.

Safety instructions before use

- Before use, check that the remote control and actuator are properly connected to the control unit.
- Check the condition of the control box and the battery.

- Check that the actuator is properly secured

- To turn on the control box, turn the emergency stop button a quarter turn to the right.
- Check that the battery is charged using the display located on the battery. Pressing and holding the BATT button displays the charge level.

Using the device

Stand assist devices allow patients with average muscle tone to be transferred to a standing position from a bed, a chair, the toilet, or an armchair in a simple and effective way.

It must be used by a caregiver with an appropriate sling.

It is important to select a sling that is suitable for the patient according to their condition and morphology. Using appropriate equipment significantly reduces the risk of injury for the patient and the caregiver.

The use of a stand assist device requires an assessment of the patient's physical capacity.

TRANSFER FROM A CHAIR


TRANSFER FROM A BED

Non-contractual visuals

Usage recommendations:

- For effective and safe standing, the patient must place their hands on the handlebar provided for this purpose.

- In addition, the patient's shins must remain in contact with the knee support at all times.
- The positioning of the knee support must be carried out by a caregiver.

Using the device

- Your stand assist device is designed to transfer a patient; it should not be used for other purposes **nor to move a patient over long distances.**
- Check that the patient's weight does not exceed the maximum weight the stand assist device can support.
- Maneuver the stand assist device by pushing or pulling the handlebar, never by pushing the patient.
- The stand assist device must be handled with care when transferring a patient and at a speed appropriate to the situation.
- Move the stand assist device only on flat and smooth surfaces. It is not recommended to use it on a slope of more than 5°; if you must use a ramp, it is advisable to have a second person assist you. Never use a stand assist device in a shower that is in operation.

The batteries of a stand assist device should never be recharged near a bath or a shower.

- Stand assist devices are useful and effective when a patient has a certain degree of dependency and needs a caregiver to perform basic care.
- They contribute to stimulating the patient and their mobility; stimulate the cardiac system; help prevent osteoporosis and disorders associated with immobility such as falls during transfers; stimulate brain activity, thus promoting the maintenance of urinary and fecal continence; improve bowel activity and bladder function.

**WARNING:
DO NOT USE
THE DEVICE IN
THE ABSENCE OF
KNEE SUPPORT**

- To use a stand assist device, it is necessary to assess the patient's physical abilities. When using a stand assist device, patients can be divided into two categories:
 1. Patients with low muscle tone
 2. Patients with muscle tone

Using the device

Braking device:

The brakes are an important safety feature of the stand assist device. The rear wheels are equipped with brakes that allow you to lock the wheel by pressing your foot on the red pedal of the wheel until it stops (wheel with brake engaged). To unlock the wheels, push the red pedal upwards (wheel with brake released).

When locking the wheels, always make sure that both wheels are locked.

Be careful when using the device on a sloped floor; both rear wheels must be locked to prevent the risk of the device moving or tipping over.

ENGLISH

COMPACT MODEL

Small wheel
front

Rear wheel
braked

ECP MECHANICAL AND ELEC- TRICAL MODEL

Front wheel
not braked


Rear wheel
braked


SUOMI

Use of the motorization

Control unit for versions
Compact and mechanical ECP


Control unit for version
Electric ECP

Remote control for versions
Compact and mechanical ECP

Remote control for version
Electric ECP

Use of the motorization

Switching on/off the control box:

- To turn on the control box, turn the emergency stop button a quarter turn to the right. The LCD screen will then display, indicating the battery charge level.
- Press the emergency stop button **after use** or **in case of imminent danger**.

Lift / lower:

- Holding down the up or down button continuously operates the arm. Releasing the button stops the movement, as does reaching an upper or lower end of travel.

Opening / closing of the legs (ECP version):

- Continuous or prolonged pressing of the open or close button operates the device's legs. Releasing the button stops the movement as well as when an open or closed end of travel is reached.

Charging the battery:

The battery block is charged directly via the control box. To do this, plug the power cord into the mains socket. Wall charger optional.


Removal / replacement of the battery:

Wall chargers are sold as optional equipment.

To remove the battery, press the gray button and pull the block upwards.

To put the battery back in place, keep the gray button pressed and position the block starting from the bottom, then release the gray button to snap the upper part into place.


Use of the sling

Nausicaa Médical recommends using the PREMIUM and VERTICONFORT back straps. All our stand assist sling straps are compatible with Nausicaa Médical stand assist devices. For advice on choosing the product best suited to your needs, please contact our sales department. To select the size of a strap, the caregiver must match the chest circumference (underarm strap) or waist (back strap) with the measurements shown in the documentation.

Size guide

SIZES	S	M	L	XL
CHEST CIRCUMFERENCE	70/90	90/115	115/140	140/165
BETWEEN THE WAIST AND THE HIPS	70/90	90/115	115/140	140/165

Before using a stand assist device, the following parameters must be taken into account depending on the situation: physical disability, the pathology and the general morphology of the patient (weight and height of the patient). Nausicaa Médical has implemented a color code to quickly identify the size of each strap:

- size S corresponds to the color Yellow
- size M corresponds to the color Red
- size L corresponds to the color Green
- size XL corresponds to the color Blue
- size TU (One Size) corresponds to the color Black

Weight / Size	1m50	1m55	1m60	1m65	1m70	1m75	1m80	1m85	1m90	+ 1m90
40 kg	S	S	S	S	S	S	M	M	M	M
45 kg	S	S	S	S	S	S	M	M	M	M
50 kg	S	S	S	S	S	S	M	M	M	M
55 kg	S	S	S	S	S	S	M	M	M	M
60 kg	M	M	M	S	S	S	M	M	M	M
65 kg	M	M	M	S	S	S	M	M	M	M
70 kg	M	M	M	M	M	M	M	M	M	M
75 kg	L	L	M	M	M	M	M	M	M	M
80 kg	L	L	L	L	M	M	M	M	M	M
85 kg	L	L	L	L	M	M	M	M	M	M
90 kg	L	L	L	L	L	L	M	M	M	M
95 kg	XL	XL	XL	XL	L	L	M	M	M	M
100 kg	XL	XL	XL	XL	L	L	L	L	L	L
110 kg	XL	XL	XL	XL	L	L	L	L	L	L
120 kg	XL	XL	XL	XL	XL	XL	L	L	L	L
130 kg	XL	XL	XL	XL	XL	XL	XL	XL	XL	XL
160 kg	2 XL	2 XL	2 XL	2 XL	2 XL	2 XL	XL	XL	XL	XL

IMPORTANT RECOMMENDATIONS

- All our strap models are Class I medical devices, they adapt and are compatible with all our stand assist devices.
- In order to ensure maximum efficiency when using the products in this range, it is essential to:
 - choose the appropriate size according to the patient
 - adjust the product as best as possible on the patient
- These products must not come into direct contact with damaged skin.
- Read the instructions for use before using the straps and carry out the checks mentioned before each use.
- Before transferring a patient with the stand assist device, make sure that the strap is properly attached to both hooks of the stand assist device.
- Sling lifespan: 3 years.

Use of the sling

Placing the strap:

- Position the bottom of the backrest at the lumbar level so that the attachment system of your strap is accessible on both sides of the patient.
 - Pass the backrest straps under the person's arms.
- Attach both straps to the hooks of the stand assist device.

Adjustments:

- The attachment system, with multiple suspension points, allows you to adjust the person's position.

Patient standing up:

- The strap is positioned at the waist, the attachment adjustment is made with the strap under tension and the patient's arms positioned at 120°.
- Hands are placed on the handles and arms are extended to activate the tone of the torso.
- The rapid and complete deployment of the lifting arm results in a standing position.

Reminder of ISO 10535:2021 standard on control:

B.2.4 Inspection of body support elements (= slings)

A periodic inspection of the flexible support element should be carried out at the frequency indicated by the manufacturer, at least twice a year.

It may be necessary to carry out more frequent inspections when a flexible support element is used or cleaned more often than usual.

Inspections should be carried out by a person with the required qualifications who is familiar with the model of the support element, as well as its use and maintenance procedures.

The inspection should allow for the detection of signs of deterioration, wear, or potential failure, and verify the legibility of the labels.

The inspection report should be kept in a secure place so it can be consulted in case of an incident.

This report should include the following information:

- the date of the inspection
- identification details of the support element and its serial number
- information regarding the condition of the support element
- the next scheduled inspection date
- the name, contact details, and signature of the inspector.

Use of the sling

Selection guide: stand assist slings:

- Stand assist sling Premium (ref.: swyp)

Recommended for patients with partial lack of tone in the legs and torso.

This sling can be used with our stand assist devices EASYLEV, WAYUP, and WAY UP GO. It can be used to assist with walking using the WAYUP. It features an under-thigh strap to allow for a semi-seated transfer.

- Stand assist sling Verticonfort (ref.: sverticonf)

Recommended for patients with partial or complete lack of tone in the legs and torso.

This sling can be used with our stand assist devices EASYLEV, WAYUP, and WAY UP GO.

Storage

To optimize storage, the Easylev Compact or ECP is designed to be folded for storage in a vertical and upright position.

To do this, it is necessary to follow the steps described below:



Step 1: Raise the lifting arm.

Step 2: Slightly unscrew the knee support and remove it from its position.

Step 3: Insert the knee support into the designated holder and screw the knob back in, then unfasten the knee support strap.

Step 4: Loosen the knobs on the upper axis of the base **WITHOUT REMOVING THEM FROM THE AXIS** which will serve as the pivot for the mast during folding.


Step 5: Completely unscrew one of the two knobs on the lower axis of the base, remove the support washer. Slightly move the mast with one hand to remove the axis with its knob and support washer with the other hand. **Always keep one hand on the mast during this step to prevent it from falling forward.**

Storage

Step 6: Tilt the upper part forward to fold the stand assist device in half, then reposition the lower axis at the end of the mast (the knobs will serve as floor supports in the upright position) with its support washers on each side of the mast and screw the knob back in place without tightening.

Step 7: Pass the knee support strap under the frame of the device to keep it folded.

- The stand assist device can now be lifted using the arm as a handle.
- It can be rolled like a trolley or a cart to take it to its storage location.
- The device is designed to be stored in a vertical position (a cover is available as an option).

Dimensions compact version
L*W*H = 580 x 570 x 1040

Dimensions ECP version
L*W*H = 600 x 520 x 1000

Reuse of the stand assist device

The stand assist device can be reused.


The stand assist device should be cleaned, disinfected, and maintained after each use.

Before any reuse, refer to the following chapters:

- Preventive Maintenance & Safety Checks
- Cleaning & Maintenance

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Cleaning & Maintenance

Objective:

- Restore the stand assist device and prevent the transmission of germs from one patient to another.
- Remove all organic soiling by mechanical action or chemical action (disinfection).

Before any handling:

- Unplug the mains power cable.
- Check that all electrical components are connected to each other.

Cleaning procedure:

- Clean the surfaces with a damp cloth using an appropriate detergent.
- Regular cleaning is recommended and should therefore be added to internal tasks.

Daily maintenance / between each patient:

- Unplug the mains power cable.
- Clean using a cloth soaked in a surface disinfectant solution.
- Clean and then disinfect all parts in direct contact with the patient.

Specific maintenance by service providers after removal of the stand assist device from the facility:

- Bio-cleaning operation.
- Steam cleaning of the various flat surfaces. Regularly change the cleaning surfaces to avoid any water load. Steam cleaning of hard-to-reach parts. For tubes, use steam with a microfiber cloth. Do not direct steam directly onto electrical boxes.

Complete cleaning and disinfection:

- Unplug the mains power cable.
- Clean the entire chassis of the device with a pre-disinfectant detergent.
- Clean the motorization with a cloth slightly dampened with a pre-disinfectant detergent solution (never clean the device, especially the electrical system, with a high-pressure cleaner, water hose, or similar.)
- Proceed with disinfection of the device (by spraying surface disinfectant or misting disinfectant solution, for example), making sure to respect the contact times of the product used.

Maintenance after cleaning and disinfection:

- Unplug the mains power cable.
- Check that the device shows no visible damage.
- Check that no parts are missing.
- Check the proper functioning of the wheels and that nothing is affecting their rolling (hair, etc.).
- Check the proper functioning of the controls as well as the connection of the remote control and actuator to the battery block.
- Clean the sockets and control buttons with a dry cloth (never clean the device, especially the electrical system, with a high-pressure cleaner, water hose, or similar).
- Check the integrity of the electrical cables (actuator and remote control).

WARNING:

- **The detergents used must have a neutral pH.**
- **Avoid abrasive products and solvents, as they could damage the surface condition of the device.**
- **Disinfection of actuators, electrical boxes, and remote controls with a microfiber wipe impregnated with disinfectant product.**

If you notice a malfunction, contact the after-sales service:

Stand Assist Devices & Patient Lifters

Phone: 04 66 71 71 80

Fax: 04 66 71 71 81 Email: sav@nausicaa-medical.com

Preventive maintenance and safety checks

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
Signature prestataire :

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1/1

Preventive Maintenance & Safety Checks

Upper part common to all 3 models


Preventive Maintenance & Safety Checks

EASYLEV Compact

ENGLISH

EASYLEV with mechanical leg spread

Preventive Maintenance & Safety Checks

ENGLISH

Document of preventive maintenance
to download from our website and to be completed.

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Signature prestataire :

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1/1

Technical & Dimensional Characteristics


Technical characteristics:

- Steel frame
- Laser-cut and bent sheets, tubes, and profiles
- Oven-cured epoxy paint
- Ball bearing wheels
- Total weight compact version: 28 kg
- Total weight ECP mechanical version: 32 kg
- Total weight ECP electric version: 35 kg
- Weight of the upper part: 16.40 kg (with one axle and two knobs)
- Weight of the base "Compact" version: 11.60 kg (with the other axle and two knobs)
- Weight of the base ECP mechanical version: 15.60 kg (with the other axle and two knobs)
- Weight of the base ECP electric version: 14.60 kg (with the other axle and two knobs)
- Maximum load: 160 kg
- Product lifespan: 8 years* (excluding electrical parts)

* Subject to regular maintenance

Dimensional characteristics:


Fixed model

Dimensions

- A. Total length: 840 mm
- B. Minimum useful height: 630 mm
- C. Maximum useful height: 1580 mm
- D. Total height: 1138 mm
- E. Frame height: 80 mm
- F. Minimum base width: 390 mm
- G. Maximum base width: 605 mm
- H. Turning diameter: 890 mm
- I. Ground clearance: 45 mm

ECP mechanical and electric model

Dimensions

- A. Total length: 1012 mm
- B. Minimum useful height: 640 mm
- C. Maximum useful height: 1550 mm
- D. Total height: 1181 mm
- E. Frame height: 115 mm
- F. Minimum base width: 570 mm
- G. Maximum base width:
 - mechanical spread: 810 mm
 - electric spread: 800 mm
- H. Turning diameter: 1120 mm
- I. Ground clearance: 75 mm

Spare Parts

**For all spare parts requests, contact
the after-sales service:**

Stand Assist & Patient Lifters

Phone: 04 66 71 71 80

Fax: 04 66 71 71 81 Email: sav@nausicaa-medical.com

Technical Specifications Motorization

Battery:

- Removable 2.2Ah Li-ion battery pack
- Digital autonomy display
- IP65

Control box:

- ABS housing
- Overload protection by thermal cut-outs
- Electronic protection in case of overload
- Emergency stop button
- Electromagnetic compatibility: compliant with EN 60601-1-2 standard
- Compliant with EN 60601-1 standards
- Protection class: II
- IP54

Actuator:

- Electric safety descent
- Low voltage DC motor 24VDC
- Power 24 V/120 VA
- Maximum push force: 3,000 N
- End-of-stroke stops by contactor
- Lifting range: 230 mm
- Lifting stroke: Strap hook / floor: Min height: 760 mm and Max height: 1690 mm. Stroke 930 mm.
- Maximum push force of the leg spread actuator: 2,000 N
- Actuator noise level: less than 60 dB at 1 meter
- IP55

Remote control:


- Remote control 2 functions / 4 functions
- Remote control with coiled cable, low voltage 5VDC
- Control button pressure: less than 1N
- IP65

Device for
class II


Wiring diagram

Technical Specifications Motorization

Control box display:

The control box has a display that allows various information to be communicated to the user:

A- Overload indicator: The load to be lifted is too heavy and may damage the device and endanger users. It is necessary to relieve the patient lift.

B- Abnormal usage duration indicator: The stand-up lift has been used for too long compared to its inactivity time. This is a safety feature intended to protect the integrity of the patient lift. The intended duty cycle is 10% operating time for 90% rest time. Let the device rest before attempting to reuse it.

C- Battery level indicator.

D- Low battery level indicator. Risk of battery damage, which must be recharged before any further use. An audible signal warns the user as soon as the battery reaches 25% of its charge.

E- Maintenance or actuator replacement indicator. The device has been used for a total duration of 600,000 seconds, or about 10,000 cycles. A general overhaul of the entire device is necessary.

Troubleshooting Guide

Symptoms	Causes	Solutions
The device, when not loaded, does not always descend.	Our devices need to have weight applied to them in order to descend.	Press down on the lifting arm while pressing the descent button on the remote control.
The moving parts of the device are stiff, the device is hard to handle.	This is due to a lack of lubrication of the moving parts.	Lubricate the moving parts.
The actuator does not work but a "click" is heard in the control box when pressing the remote control.	1. The batteries are discharged. 2. The actuator cable is not connected. 3. The actuator cable is damaged. 4. The actuator or the control box is faulty.	1. Charge the device. 2. Plug in the actuator cable. 3. Proceed with replacing the actuator. 4. Have the entire electrical system serviced.
The actuator does not work but no "click" is heard in the control box when pressing the remote control.	1. The emergency stop button is pressed. 2. The battery is discharged. 3. The remote control cable is not connected or is poorly connected. 4. The remote control cable is damaged. 5. The entire electrical system is damaged.	1. Unlock the emergency stop button. 2. Check the condition of the charger and then charge the battery. 3. Plug in the remote control. 4. Replace the remote control. 5. Have the entire electrical system serviced.
The actuator stops.	1. The battery is discharged (a beep sounds) 2. The device has been used for too long. 3. Actuator poorly connected	1. Check the condition of the charger and then charge the battery. 2. Let the device rest for a while (protection for the actuator's lifespan – 2 min of use for 18 min of rest). 3. Connect the actuator
The actuator does not go up or down anymore.	1. The emergency stop button is pressed. 2. The box is faulty. 3. The remote control is faulty. 4. The battery is empty	1. Unlock the emergency stop button. 2. Have the box serviced. 3. Replace the remote control. 4. Recharge the battery

Warranty

• Article 1: NAUSICAA Médical S.A.S. guarantees this device against all manufacturing and assembly defects of its mechanical and electrical components, and this only for devices used under the conditions specified by NAUSICAA Médical S.A.S.

The warranty covers the mechanical and electrical parts.

This warranty, the conditions of which are defined below, is valid for 5 years except for the batteries (lithium battery: 2 years warranty, lead battery: 1 year warranty).

• Article 2: The warranty entitles you to free labor as well as the free replacement of parts found to be defective.

• Article 3: The outbound shipping of the device, as well as all related costs, are the responsibility of the reseller. The goods always travel at the risk and under the responsibility of the reseller.

Under warranty: the return costs after intervention will be borne by NAUSICAA Médical S.A.S.

Out of warranty: return costs are at the reseller's expense whether or not they accept the repair estimate.

• Article 4: The warranty does not apply if the claims result from:

- an accident, misuse of the device, or negligence by the buyer.
- transportation of the device carried out without adequate protection.
- a modification or transformation not validated by NAUSICAA Médical S.A.S.
- the impact of external agents (natural disaster, fire, shocks, humidity, flooding, lightning, etc.).
- installation and/or use in a manner not compliant with technical and safety standards in the event that the device operates in a country other than the country of purchase; and/or if the power supply is not suitable for the operating voltage of the device.
- a lack of regular maintenance.

• Article 5: The reseller cannot claim the benefit of the warranty:

- if the serial number of the device has been removed, modified, or made illegible.
- if the device under warranty has been modified without the approval of NAUSICAA Médical S.A.S.

• Article 6: During the repair of defective equipment, no equipment will be loaned.

• Article 7: Any warranty claim must be made through the reseller. If this is not possible, the buyer may possibly send their equipment directly to NAUSICAA Médical S.A.S. In this case, the buyer must include with the equipment a letter stating the reseller's contact details and a copy of the purchase invoice.

• Article 8: The sending of spare parts under warranty will only be done after consultation with the After-Sales Service of NAUSICAA Médical S.A.S.

Please note that defective spare parts must be returned to the After-Sales Service of NAUSICAA Médical S.A. or they will be invoiced 1 month after the parts are sent by NAUSICAA Médical S.A.S.

• Article 9: Defective parts replaced under or out of warranty will themselves be guaranteed for 6 months from the date of repair or the sending of spare parts.

• Article 10: No reseller may unilaterally modify the terms of this warranty.

NAUSICAA

MEDICAL

Head Office Sales Department
Phone: 04 66 51 50 80
Fax: 04 66 51 50 47
Email: contact@nausicaa-medical.com
www.nausicaa-medical.com

Your contact:

Manufactured by NAUSICAA Medical S.A.S.

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- muutos tai muunnos, jota NAUSICAA Médical S.A.S. ei ole hyväksynyt.
- jos takuulla olevaa laitetta on muutettu ilman NAUSICAA Médical S.A.S.:n hyväksyntää.
Valmistaja: NAUSICAA Médical S.A.S.