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V03.09/26

This device is intended for mobile patient lifts allowing to transfer a person in standing position.

Before each use, it is important to check the presence of the label, the stitches of the slings, that they do not fray or break, that the device does not show any tears, snags, holes or signs of fraying (fabric and slings), that the buckles are not cracked or have no weak fasteners

PADDED BUTTOCKS COMFORT BELT


INSTALLATION

Clip the 50 mm male-female buckles to one side of the sling, pass the comfort belt under the patient's buttocks, clip the male-female buckle to the other end of the sling, check that the patient is correctly installed in the under the buttocks sling before doing a transfer.

Extract from the HAS national commission for the evaluation of medical devices dated 14/09/2010: Slings are accessories that must be changed regularly (hygiene, bathing, etc.). Two slings per patient must be covered on purchase and annually if necessary. Each sling must be dedicated to a single patient.

The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative. This medical device must be implemented or trained by a health professional.

SFC-PF

V03.09/26