SANGLES

NAUSICAA MEDICAL

# **Ceinture** mini pelvienne

## Descriptif

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## Description

Pelvic belt made from high-strength polyester mesh. Prevents the pelvis from sliding forwards.

## Recommendations

The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.

## Composition

• 100% polyester mesh
• 100% polyester mesh
• 100% polyester bias binding
• 100% polyester straps
• Nylon buckles

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CE

## Sizes and dimensions

It is made from high-strength polyester mesh.
• Available in 3 sizes: S, M, L
• Maximum usable circumference: 230 cm
• Waist circumference (in cm)

## Fitting

This medical device is intended to ensure patient safety.
Before each use, it is important to check: that the label is present; the stitching on the straps: ensure they are not fraying or coming apart; that the device has no tears, snags, holes or signs of fraying (fabric and fastening straps); and that the buckles are not cracked or showing any signs of weak fastenings.

• Position the pelvic belt section on the seat of the chair, with the narrow end against the backrest.
• Place the narrow strap around the backrest, adjust it firmly and fasten the quick-release buckle (see front).
• Have the patient sit in the chair.
• Pass the pelvic belt between the patient's legs and pull it up onto the abdomen.
• Adjust the pelvic section to avoid any creases.
• Pass the remaining strap behind the backrest, fasten the buckle and adjust the length, taking care not to compress the abdomen.

## Important recommendations

• To ensure maximum effectiveness when using the products in this range, it is essential to:
- choose the correct size for the patient
- fit the product as closely as possible to the patient
• These products must not come into direct contact with broken skin.

## End-of-life management

• Medical devices at the end of their useful life (obsolete, beyond repair or past their service life) must not be disposed of with household waste.
• The user is responsible for decommissioning them, cleaning and disinfecting them, sorting the materials and directing them to the appropriate treatment channels, in accordance with current regulations.
• Textile components (polyester, polyamide, polyurethane, nylon) and plastics must be recycled where possible.
• Used or damaged slings must no longer be used for lifting.
• Contaminated devices must be disinfected and then disposed of in accordance with applicable procedures, in particular the DASRI regulations for infectious waste.
• Material recovery and recycling must be prioritised to ensure the sustainable, safe and compliant management of medical devices at the end of their life.