Ceinture pelvienne confort

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Descriptif

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FOR ALL CHAIRS

Description

Pelvic belt with internal steirthane protection waterproof and breathable.

- Available in 3 sizes S, M, L
- Maximum useful area: 230 cm

Sizes

Waist size (in cm)	
S	50/80
M	70/100
L	90/130

CPCF3HP-PF

This medical device must be handled only by employees trained (by a healthcare professional) and qualified for its use (intended purpose, installation and environment of use, etc.). If the medical device is incorrectly installed, there is a risk of accident for the patient.

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Recommendation

The use of the device requires the presence and vigilance of a third person. This medical device must be implemented or trained by a health professional.

The comfort pelvic belt prevents the pelvis from sliding forward. This product is intended exclusively for the medical use specified above.

Composition

- Combination of polyurethan foam and on both faces jersey/PU coating.
- 100% polyester bias
- 100% polyester strap
- Nylon buckle

Set up

This medical device is intended for patient holding and positioning. It is not a restraint system.
Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

- Place the part « pelvic holder » on the seat of the chair, the narrowest part to the backrest.
- Pass the strap of the narrowest part around the backrest of the chair, firmly adjust it and close the buckle.
- Sit the patient on the chair and raise the belt part between his legs and adjust the belt around the patient's abdomen.
- Adjust the pelvic part to avoid creases.
- Pass the strap which has not been used yet behind the backrest, close the buckle, adjust the length of the strap while being careful not to compress the abdomen of the patient.

Important recommendations

- In order to get a maximum efficiency of these sort of product, it is necessary :
- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient
- These products must not be in direct contact with a wounded skin.

End-of-life management