NAUSICAA
MEDICAL

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Made in France
since 1992

www.nausicaa-medical.com

This medical device is intended for patient positioning.This is not a restraint system.
Before each use, it is important to check: the presence of the label, that the device does not show any tears, snags, holes or signs of fraying and the good condition of the hook and loop.

FOOT SUPPORT

TO USE ONLY FOR THE INDICATED USE :
KEEPING YOUR FOOT ON THE WHEELCHAIR PALLET

The foot support keeps the foot on the pallet of the chair (with self-gripping attachment under and on the sides of the pallet).

⚠

Restraint is exceptional and must be limited in time (requires a medical prescription). It must be accompanied by a regular assessment of the benefits/risks and increased patient monitoring. It should only be used by staff who are trained and qualified in its use. In the event of incorrect use, there is a risk of accident to the patient.

⚠

The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative. This medical device must be implemented or trained by a health professional.

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CE

Important Recommandations :

• In order to get a maximum efficiency of these sort of product, it is necessary :

- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient

• These products must not be in direct contact with a wounded skin.

FOOT SUPPORT

HOSPITAL RANGE
Polyurethane-coated jersey lined with polyurethane foam
Bias 100% polyester
Self-gripping 100% polyester
PU foam 100 kg / m3

• Attach each foot support to each wheelchair paddle with the self-gripping sole of the foot support.

• Place each foot of the patient in a footrest.

• Close the footrest on the patient's foot with the self-gripping device.

End-of-life management