NAUSICAA
MEDICAL

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Made in France
since 1992

www.nausicaa-medical.com

This medical device is intended for patient support and positioning. It is not a restraint system.
Before each use, it is important to check the presence of the label, the stitches of the straps, that they do not fray or break, that the device does not show any tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

SECURITY MITTS

USES :
COMFORTABLE, THEY PREVENT THE RESTLESS PATIENT FROM TEARING HIS CARE DEVICES AS WELL AS HIS INCONTINENCE PROTECTIONS.

Security mitts prevent the restless patient from tearing his care devices as well as his incontinence protections.
The breathable net part on the top increases comfort and the integrated aluminium stiffener prevents the hand flexion.

⚠

Restraint is exceptional and must be limited in time (requires a medical prescription). It must be accompanied by a regular assessment of the benefits/risks and increased patient monitoring. It should only be used by staff who are trained and qualified in its use. In the event of incorrect use, there is a risk of accident to the patient.

⚠

The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.
This medical device must be implemented or trained by a health professional.

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CE

Important Recommandations :
• In order to get a maximum efficiency of these sort of product, it is necessary :
- to get the best possible adjustment to the patient
• These products must not be in direct contact with a wounded skin.

SECURITY MITTS

HOSPITAL RANGE
Combination of polyurethan foam
and on both faces jersey/PU coating
100% polyester net
100% polyester bias
100% polyester strap
Nylon buckle
Aluminium stiffener

- Place the patient's hand in one of the security mitt, net part on the upper part of the hand.
- Close the nylon buckle, then adjust the tighthening.
- Secure it with the velcro strap.

End-of-life management