Sangle de rehaussement

Descriptif

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Raising sling

Description

This sling allows the caregiver to raise the bed alone, even with uncooperative patients. It is easy to install and protects the caregiver's back.

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Recommendation

Use of the device requires the presence and vigilance of a third party.
This medical device must be operated or trained by a healthcare professional.

This technical aid is designed for caregivers to help them lift a patient with a low level of physical activity into bed with less effort and without damaging their backs.

Composition

- Sterithane laminate
- 40 mm polyester webbing – red
- 22 mm Polyester bias - black
- Polyamide handle

Recommendations for use

Before each use, it's important to check: that the label is attached; that strap stitching is not loose or frayed; that the device is free of tears, cracks, holes or signs of fraying (fabric and hooking straps); and that buckles are not cracked or showing signs of weak attachment.

Single-person use:
The raising sling should be hooked to one side of the bed at a fixed point towards the top of the bed (slatted frame, gallows). The padded sleeve should be positioned at the patient's buttocks. The caregiver retrieves the other end of the sling from the other side of the bed and uses a splitting movement to raise the patient. If able, the patient can participate by looking at his stomach and pulling on the gallows handle and/or pushing on his feet.
Two-person use:
The sheath is placed on the patient's buttocks, with one person on each side of the bed holding the strap by its handles. The two people synchronize their movements to lift the patient to the top of the bed. As mentioned in the previous paragraph, the patient can participate to the extent of his or her ability.
Use with sliding sheet :
The raising sling can be used in combination with a sliding sheet, which is fitted beforehand. This combination enables effortless raising of heavier patients or those with less physical capacity to participate.

End-of-life management

- Medical devices at the end of their useful life (obsolete, beyond repair or past their service life) must not be disposed of with household waste.
- The user is responsible for decommissioning them, cleaning and disinfecting them, sorting the materials and directing them to the appropriate treatment channels, in accordance with current regulations.
- Textile components (polyester, polyamide, polyurethane, nylon) and plastics must be recycled where possible.
- Used or damaged slings must no longer be used for lifting.
- Contaminated devices must be disinfected and then disposed of in accordance with applicable procedures, in particular the DASRI regulations for infectious waste.
- Material recovery and recycling must be prioritised to ensure the sustainable, safe and compliant management of medical devices at the end of their life.