NAUSICAA
MEDICAL

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Made in France
since 1992

www.nausicaa-medical.com

This medical device is intended for patient positioning.This is not a restraint system.
Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

HEADREST

The removable partial head support device is used with the following harnesses:SURPC2-PF, SUP2-PF, SURP2-PF, HA2-RBA-PF, SUBAS2-PF, SUBAI-PF, SUCBAS2-PF, SUCBAI-PF, SST-PF, SSTBA-PF, SHA2-PF, SHAC2-PF, SHA2BA-PF, SHAC2BA-PF, SUCPEHA2-PF, SUCBAHA2-PF.
Slide the tip upward between the patient's back and the strap. Once the strap is tightened, the headrest will be secured.

TAMOV-PF ☐

The removable, ribbed partial head support device is used with the following slings: STORP2-PF, SATL-PF, STO-PF.
Slide the narrowest part downward between the sling and the patient. Once the sling is tightened, the headrest will remain in place.

TABAL-PF ☐

The TAAT-PF headrest is used with the following harnesses: SUBAS2-PF, SUBAI-PF, SUCBAS2-PF, SUCBAI-PF, SHA2-PF, SHAC2-PF, SHA2BA-PF, SHAC2BA-PF, SUCPEHA2-PF, SUCBAHA2-PF
It connects to the straps via the two 50-mm strap ends located on either side of the handle at the top of the strap near the nape of the neck. Once secured, the headrest becomes an integral part of the strap and attaches to the same hooks as the shoulder straps.

TAAT-PF ☐

⚠

Remove the aluminium bar before putting the headrest in the washing machine.

TAMOV-PF / TABAL-PF / TAMAX-PF

⚠

The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.
This medical device must be implemented or trained by a health professional.

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