



PROTECTIVE COVER



This cover, to be sterilised before use, protects the lift after it has been disinfected. It has a front pocket to insert documents. To optimise storage, the lift can be folded and stored in an upright position with a protective cover (refer to the 'Storage' section of the user manual for step-by-step instructions on folding the lift).


NB: It is possible to roll the lift like a trolley or caddy to place it to its storage location.

Available for the stand up models below :

- Way Up'Go Compact and Way Up'Go ECP
- Blue Wayup 3 Compact and Blue Wayup 3 ECP

Composition
75% Polyamide
25% Polyurethane foam

Dimension
Height : 95 cm
Width: 62 cm
Deep : 57 cm



HPD-WP

The stand up lift and sling should be cleaned, disinfected and maintained after each use, before storage and before re-use following storage.



v 02.08/26