

# SELECTION GUIDE: SLINGS

Nausicaa Médical slings can be used with any lift or transfer rail fitted with a smooth hook system, open or closed, capable of receiving a 40 mm wide sling without impacting the strength of the sling on beams ranging from 2 to 4 points, which excludes hooking with a metal plate with a hole. In accordance with the ANSM document “Compatibility between medical devices” of July 2007: Nausicaa Médical sit-to-stand slings comply with standard EN ISO 10535: 2021 and our risk management process demonstrates that the use of Nausicaa Médical slings with sit-to-stand slings fitted with a 2-point hooking system does not give rise to any additional risks compared with the original device.

	RECOMMENDED SLINGS FOR TRANSFERS		Slings recommended for toileting. Cannot be left under the patient.
	HAMMOCK SLINGS (Can remain under the patient)	U Shape Slings (Remove after use)	
Loss of tone in the lower limbs with chest support and with head support	Hammock sling	U shape sling U shape sling with quick adjustment	Toilet sling Quick transfer U shape sling
Loss of tone in the lower limbs without chest support and with head support	Hammock sling	U shape sling U shape sling with quick adjustment	Toilet sling
Loss of tone in the lower limbs without torso support and without head support	Hammock sling with head support	U-shaped backrest sling Quick-adjust U-shaped sling with headrest option Ergonomic U-shaped sling	
Amputee patient	Hammock sling without crossing		
Agitated/anxious patient	Hammock sling without crossing	Wrap around U shape sling U shape sling with backrest	
Spastic patient	Hammock sling without crossing and with head support	U shape backrest sling Quick-adjust U shape sling with headrest option Ergonomic U shape sling	Toilet sling with headrest option U shape sling with quick adjustment
Physically and mentally dependent patient	Hammock sling	Wrap-around U shape sling U shape sling with backrest Quick-adjust U shape sling with headrest option Ergonomic U shape sling	
Contraindication for hip constriction	Hammock sling without crossing	Fast transfer sling	Fast transfer sling



All our textiles are REACH compliant.