

# PELVIC SUPPORT BELT



CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

The bed support belt is designed to prevent falls and provide safe positioning for the patient. It is compatible with the NAUSI-SURE magnetic closure system, for beds up to 90 cm wide.



REFERENCES	
CMLP2HP-PF-	
S	ou M ou L ou XL

Wrist size (cm)

S	60 / 90
M	80 / 105
L	95 / 125
XL	115 / 145

# FASTENERS



Wrist straps



Ankles straps

CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

Wrist straps are fitted around the patient's wrists to prevent bruising and injury, for example from infusions being pulled out. A safety version is available for restraining highly agitated patients (ref.: LMATPOS2-PF). Ankle straps keep the legs slightly apart to avoid the risk of chafing.



REFERENCES	
LMATPO2-PF (2 units)	
LMATPOS2-PF (2 units)	
LMATCH2-PF	