


SLIDE SHEETS



CONFORT STYLE BIEN ÊTRE
FABRICATION FRANÇAISE

For patient handling and raising, the sliding sheet is a device that makes it easier to move a patient lying in bed, reducing the risk of back pain and musculoskeletal disorders for carers. Available in two widths: 70 and 110 cm. Available in different lengths: width 70: 80, 130 and 195 cm / width 110: 140 cm / width 48: 48 cm.

SLIDE SHEETS WITH HANDLES



CONFORT STYLE BIEN ÊTRE
FABRICATION FRANÇAISE

For patient handling and raising, the sliding sheet is a device that makes it easier to move a patient lying in bed, reducing the risk of back pain and musculoskeletal disorders for carers. The handles make it easier to grip. Available in 2 widths: 70 and 110 cm. Available in several lengths: width 70: 80, 130 and 195 cm / width 110: 140 cm

REF. 1 PIECE	REF. 1 PIECE	REF. 5 PIECES	REF. 5 PIECES
DGL-195-70-PF DGL-130-70-PF DGL-80-70-PF	DGL-140-110-PF DGL-48-48-PF	DGL5-195-70-PF DGL5-130-70-PF DGL5-80-70-PF	DGL5-48-48-PF

REF. 1 PIECE	REF. 5 PIECES
DTR-195-70-PF DTR-140-110-PF DTR-130-70-PF	DTR5-195-70-PF DTR5-140-110-PF DTR5-130-70-PF

