

NAUSIFLOW 2S PUMP



FEATURES

- 
Manual patient comfort adjustment
- 
Seated position
- 
Low-pressure static mode for palliative care treatments in particular
- 
Low pressure indicator + warning beep. Normal pressure indicator.

The NAUSIFLOW 2S pump features a pressure adjustment knob for setting the desired comfort level. Its high-flow mechanical pump maintains a constant pressure at the chosen comfort setting.

Two operating modes:

- Dynamic mode (alternating 1 on 2)
- Palliative care mode (static low pressure)

Can be used with our various mattresses.

An audible alarm can be deactivated via the button on the front panel.

For added convenience, the air filter and fuses are easily accessible from the rear of the pump.


PREVENTIVE MAINTENANCE TRAINING
 Nausicaa Médical trains you in the care and maintenance of our products to enhance their durability and help you get the most out of your investment.

REFERENCE
 NA2S-COMP-PM

 130 kg



WARRANTY
2
 YEARS

