

ECO STAND UP SLING

- Only available in jersey
- One size fits all



The sling is positioned at waist level on the WAYUP and the adjustment is made with the sling with the patient's arms positioned at 120°.

The sling is positioned at waist level and/or sub-axillary on the EASYLEV according to the patient's tone.

HOOK ADJUSTMENT



Two levels of hooks

SLINGS POSITIONNING



UNCROSSED

Decent bust tonicity



CROSSED

Lack of bust tonicity

REF. SXPECO2-PF-TU

