

# VERTICONFORT STAND UP SLING

The sling is positioned at waist level on the WAYUP and the adjustment is made with the sling with the patient's arms positioned at 120°.

The sling is positioned at waist level and/or under the arms on the EASYLEV, depending on the patient's tone. Possibility of semi-seated transfer thanks to its sub-gluteal sling and securing the patient in walking aid mode only on the WAYUP.



- Only available in comfort jersey
- Sizes S, M, L, XL.

## HOOK ADJUSTMENT



Four levels of hooks.

## SLINGS POSITIONNING



Ref.: SWYP3-PF

With bottom sling



Ref.: SWYP3-PF

Walking aids



Ref.: SVERTICONF-P



Ref.: SWYP3-PF

CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE



REF.	SVERTICONF-PF SWYP3-PF
OPTION	SFC (Only for the SWYP model)