

SHAPE ECO SLING

- Available in jersey, comfort jersey and mesh
- One size



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE

REF.	SURPC3-PF-TU (jersey comfort) SUP3-PF-TU (jersey) SURP3-PF-TU (mesh)
OPTION	TAMOV-PF (removable headrest)



The U-shape prevents the slings from being placed under the patient's buttocks, making fitting and removal more comfortable for the patient and easier for the carer.


This sling is available in comfort jersey with comfort foam in the leggings to reduce shearing, or in mesh for bathing. This sling must not be left under the patient.

HOOK ADJUSTMENT



This Eco sling offers a single adjustment.

SLINGS POSITIONNING



HEADREST (OPTION)

The addition of the removable headrest (TAMOV-PF) helps to protect the cervicals when picking up on the ground, for example.

