

# ENVELOPING SLING

For agitated and anxious patients.

- Only available in comfort jersey
- Sizes S, M, L, XL.



CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

REF.	SLK-PF
------	--------



The integrated whalebone headrest and the enveloping shape of this sling are recommended for transfers of agitated and anxious patients. Its U-shape prevents the sling from being placed under the patient's buttocks, making application and removal more comfortable for the patient and easier for the carer.

The positioning of the hooks at shoulder level, running along the entire back of the sling, prevents compression of the patient during the transfer and is suitable for 2-point slings.

## HOOK ADJUSTMENT



The multiple hooks system allows the patient's position to be varied. For optimum comfort in the seated position, we recommend adjusting the sling as short as possible on the shoulders (red setting) and as long as possible on the legs (black setting).

## SLINGS POSITIONNING



CROSSED BETWEEN LEGS

Recommended choice



CROSSED UNDER LEGS

Second option