

# HAMMOCK ECO SLING

- Available in jersey and mesh
- One size fits all



CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

<b>REF.</b>	HA3-RBA-PF-TU (mesh) HA-PF-TU (jersey)
<b>OPTION</b>	TAMOV-PF (removable headrest)



The simple cut and scalloped sides make it the sling of choice for leaving under the patient in the armchair, as it has no excess thickness or folds which encourage pressure sores and discomfort.

It should be applied to the bed in a supine position at 180°, with the patient in a lateralization.  
The sling can remain under the patient.

## HOOK ADJUSTMENT



This Eco sling offers a single adjustment.

## SLINGS POSITIONNING



## HEADREST (OPTION)

The addition of the removable headrest (TAMOV-PF) protects the cervicals when picking up on the ground, for example.

