

# PREMIUM JERSEY HAMMOCK SLING

The simple cut and scalloped sides make this sling ideal for leaving under the patient in the chair, as there are no excess thicknesses or creases that encourage pressure sores and discomfort.

This sling is available in comfort jersey with comfort foam in the leg pads to reduce shearing.

It should be fitted in bed in a supine position at 180°, with the patient positioned laterally.

- Available in jersey and comfort jersey
- Sizes S, M, L, XL.

## HOOK ADJUSTMENT



The multiple adjustment system allows the patient's position to be varied. For optimum comfort in the semi-lying position, we recommend adjusting the sling in the intermediate position to a yellow setting for the shoulders and black (or red) for the legs.

## SLINGS POSITIONNING



CROSSED

Recommended choice



CROSSED OVER LEGS

Second option

## HEADREST (OPTION)

The headrest with fasteners (TAAT-PF) and removable headrest (TAMOV-PF) are available as options.



CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE



REF.	SHA2-PF- (jersey) SHAC2-PF- (jersey comfort)
OPTION	TAAT-PF (headrest with fasteners) TAMOV-PF (removable headrest)