ERGONOMIC U SHAPE SLING

The back section is extended by the rigid foam headrest, which is lined with 3D fabric, making this sling particularly comfortable and helping to support the spine during transfers and picking up from the floor.

The U-shaped design prevents the sling from being placed under the patient's buttocks, making it more comfortable for the patient and easier for the caregiver. It can only be fitted in the seated position.

- Only available in comfort jersey
- Sizes S, M, L, XL.

HOOK ADJUSTMENT



The large loop means you can adjust the inclination without having to unhook the sling. The multi-hanger system allows the patient's position to be varied.

For optimum comfort in the seated position, we recommend adjusting the red strap and the legs as far as possible (black adjustment).

SLINGS POSITIONNING



Recommended choice



Second option





