

MATTRESS NAUSIFLOW 834 XXL

Recommended for stage 3 and 4 pressure sores



Attention: In dynamic mode, no discharge in the sacral zone.

The NAUSIFLOW 834 is an all-air mattress with a total height of 25 cm. It is made up of 20 independent, disconnectable cell-on-cell cells, i.e. 15 cm of active air + 10 cm of confined air, ensuring relative safety in the event of damage to the upper cell. Its 4 cells, located in the heel area, are fitted with a quick connector enabling localized discharge while leaving the mattress operational by fitting a plug.

The mattress is covered by a fully-zipped cover made of liquid-proof, breathable, PU-coated bielastic jersey, treated with silver ions (Ag+) for antibacterial and antifungal properties. It is M1 fire-resistant and machine-washable up to 95°C.

Available in 120 and 140 cm width and can only be used with the 100 MAXI pump.



Nausicaa Médical trains you in and maintenance of our and maintenance of our products durability and enable you to maximize your return on investments.



REFERENCES

NA1034NPU-MAT120
NA1034NPU-MAT140