

NAUSIFLOW 512 MATTRESS

Recommended for stage 1 and 2 pressure sores



CAUTION: In dynamic mode, no discharge in the sacral zone



The NAUSIFLOW 512 mattress has a total height of 17 cm. It comprises 17 independent, disconnectable cells with an air height of 13 cm and a 4 cm thick stabilizing foam base. The 4 cells located in the heel area feature a quick connector for localized discharge, while leaving the mattress in place. The mattress remains operational when a plug is fitted. The mattress is covered by a fully-zipped cover made of liquid-repellent, breathable PU-coated bielastic jersey, treated with silver ions (Ag+) for antibacterial and antifungal properties. It is M1 fire-retardant and machine washable up to 95C°. Available in widths of 70, 88 and 120 cm, and compatible with 2S, AUTO QUATTRO and MAXI pump for 120 cm widths.



3 valves air supply.



PREVENTIVE MAINTENANCE TRAINING

Nausicaa Médical trains you in the care and maintenance of our products to improve their durability and enable you to maximize your return on investment.



The first step in preventing bed sores is to change the body's support points every three hours on average, in order to limit prolonged compression of the vessels opposite bone protrusions, and thus enable better irrigation of the various tissues. Nausicaa Médical offers a wide range of therapeutic supports, the result of years of experience in the field and listening to healthcare professionals. The aim of our range is to provide a lying surface suitable for people who are bedridden for long periods, and who have little or no mobility. Our devices can be adapted to most bed widths, and are selected according to the patient's condition, morphology, weight and the number of hours spent in bed.

REFERENCES	REFERENCES
NA512NPU-MAT70 NA512NPU-MAT88	NA512NPU-MAT120



150 kg

WARRANTY 2 YEARS

