

ECO TOILET SLING

- Only available in comfort jersey
- Universal size



The enlarged U-shape of this sling frees up the patient's buttocks completely, making it easier to position the patient on the toilet, and makes the sling easier to put on and take off.

The legs pass from the outside to the inside. This sling is positioned with the arms on the outside.


This type of sling is not recommended for patients with very low muscle tone.

CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE

REF.	STO-RP3-TU (jersey)
OPTION	TABAL-PF (removable headrest)




HOOK ADJUSTMENT



This Eco sling offers only one adjustment.

SLINGS POSITIONNING



HEADREST (OPTION)

The addition of the removable, whalebone headrest (TABAL-PF) protects the cervicals when picking up on the ground.

