

# CHEST SUPPORT

## EVOLUTIVE ECO BUST BELT



The chest moves forward.



The loops on our belts are freely positionable, so they are accessible to the patient.

CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

By assembling the bust support module on the Eco evolutive belt, you'll obtain a bust support.

- Universal size
- Maximum useful perimeter: 228 cm

## Y SLIM BUST BELT



The loops on our belts are freely positionable, so they are accessible to the patient.

CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE

The Y Slim bust belt allows separate adjustment of shoulder heights to help position and compensate for poor posture.

- Sizes S, M, L, XL
- Maximum useful perimeter: 230 cm

Waist size (cm)

REFERENCE

CSY2HP-PF

S	60 / 80
M	75 / 95
L	90 / 110
XL	110 / 130

