

PREMIUM JERSEY U SHAPE SLING

The U-shaped design prevents the slings from being placed under the patient's buttocks, making fitting and removal more comfortable for the patient and easier for the carer. The positioning of the hooks at shoulder level, covering the entire back of the sling, prevents the patient from being compressed during the transfer.

The comfort foam in the leg slings reduces shearing. This sling must not be left under the patient.

- Available in jersey and comfort jersey
- Tailles S, M, L, XL.

HOOK ADJUSTMENT



The multiple suspension system allows the patient's position to be varied. For optimum comfort in the seated position, we recommend adjusting the sling as short as possible on the shoulders (yellow setting) and as long as possible on the legs (black setting).

SLINGS POSITIONNING



CROSSED BETWEEN

Recommended choice



CROSSED BELOW LEGS

Second option

HEADREST (OPTION)

The headrest with fasteners (TAAT-PF) and removable headrest (TAMOV-PF) are available as options.



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE



REF.	SUBAS2-PF- (jersey) SUCBAS2-PF- (jersey confort)
OPTION	TAAT-PF (headrest with fasteners) TAMOV-PF (removable headrest)