

## HEEL PADS



Small Model



Large Model

The heels pads wrap the foot in a comfortable lining, preventing rubbing on the heels and malleoli.

They are available in two different sizes:

- Small (LM2TAL1)
- Large (LM2TAL2)



### REFERENCES

LM2TAL1-PF (the pair)  
LM2TAL2-PF (the pair)

## BED BARRIER PROTECTIONS



CONFORT STYLE BIEN ÊTRE  
FABRICATION FRANÇAISE

The foam protectors fit over the bed rails, preventing the patient from bumping into the rails and passing their limbs through them, and protecting them from contact with the metal. Quick and easy attachment to the bed rails using quick-release buckles. The removable protectalaise cover is M1 fire-rated.

Available in two lengths: 140 cm (small model) and 190 cm (large model). Sold in pairs only.



### REFERENCES

PROBAR-M-PF  
PROBAR-L-PF