

# BLUE WAYUP 3


## COMPACT

Its fixed base allows it **to be adapted to all spaces, even confined ones such as toilets.**



### Dimensions

- A. Total length: 91,5 cm
- B. Minimum useful height: 77 cm
- C. Maximum useful height: 116 cm
- D. Total height: 94,5 cm
- E. Frame height: 8 cm
- F. Minimum base width: 39 cm
- G. Maximum base width: 60,5 cm
- H. Turning diameter: 89 cm



### REFERENCES

WAYUPBLUE3-CBLI

# BLUE WAYUP 3

## MECHANICAL OR ELECTRICAL OPENING LEGS

Its wide mechanical or electric base legs enable it **to be used with all types of wheelchairs.**


Its removable platform allows it **to be used as a walking aid or as a lowered platform.**

**Improved rolling quality** thanks to its 100 diameter castors with double ball bearings.



### Dimensions

- A. Total length: 100 cm
- B. Minimum useful height: 74,5 cm
- C. Maximum useful height: 117 cm
- D. Total height: 99 cm
- E. Frame height: 11,5 cm
- F. Minimum base width: 57 cm
- G. Maximum base width mechanical: 81 cm  
Maximum base width electrical: 80 cm
- H. Turning diameter: 112 cm



### REFERENCES

WAYUPBLUECP3-CBLI  
WAYBLUE+3-CBLI