



This medical device is intended for patient positioning. This is not a restraint system.

Before each use, it is important to check: the presence of the label, the stitches of the straps: that they do not fray or break, that the device does not show tears, snags, holes or signs of fraying (fabric and straps), that the buckles are not cracked or have no weak fastening.

LIFTING BELT

TO BE USED ONLY FOR THE
STATED PURPOSE :
AID TO LIFT A PERSON
SITTING

Jersey/PU coating
PE foam 100 kg/m³
Polyester strap
100% polyester fabric




Description

The Lifting Belt is used to raise a sitting person and preserves the caregiver of awkward postures during manual transfers.

Its inside fabric is in anti-slipping PU coating.



CREL-PF


The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.
This medical device must be implemented or trained by a health professional.



v 05.02/23