




www.nausicaa-medical.com

Before each use, it is important to check: the presence of the label, that the device does not show any tears, snags, holes or signs of fraying.

SLIDE SHEETS



DGL-PF



The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.

This medical device must be implemented or trained by a health professional.



v 05.02/23



SLIDE SHEET

HOSPITAL RANGE

Tissu 100% polyamide :

- Enduction silicone (dessus)
- Enduction polyuréthane (dessous)

DESCRIPTION

- For the handling of a patient, the Slide Sheet is a device that facilitates lateral and vertical translation as well as the raising of a lying patient in a bed.
- The sliding effect facilitates the work of carer, promotes the prevention of pressure ulcers and improves patient comfort.
- The Slide Sheet is composed a of synthetic fabric polyamide tube in polyurethane induction, the two layers of fabric in contact slide on one another.
- Silicone Treatment on Patient Side.

USES

