

EASYLEV 5

COMPACT

Its very low ground clearance (6.5 cm) allows it to fit under very **low chairs, especially 6-wheel electric chairs.**
Its fixed base allows it **to be adapted to all spaces, even confined ones such as toilets.**




Dimensions

- A. Total length: 92 cm
- B. Minimum useful height: : 72 cm
- C. Maximum useful height: 150 cm
- D. Total height: 124 cm
- E. Frame height: 6,5 cm
- F. Minimum base width: 34 cm
- G. Maximum base width: 56,5 cm
- H. Turning diameter: 90 cm

REFERENCES

EASYLEV5-PF
EASYLEV5-LI-PF



EASYLEV 5

MECHANICAL OR ELECTRICAL OPENING LEGS

Its wide mechanical or electrical opening legs makes it **compatible with all types of chairs.**
A lowered platform is available as standard.



Dimensions

- A. Total length: 107 cm
- B. Minimum useful height: 69 cm
- C. Maximum useful height: 152 cm
- D. Total height: 130 cm
- E. Frame height: 11,5 cm
- F. Minimum base width: 54 cm
- G. Maximum base width mechanical: 83 cm
Maximum base width electrical: 80 cm
- H. Turning diameter: 110 cm

REFERENCES

EASYLEV5ECP-M-PF
EASYLEV5ECP-M-LI-PF
ESL5ECP-E2-PF / ESL5ECP-E2-LI-PF

