

PRODUCT DATASHEET

BLUE EASYLEV 3 STAND UP LIFT COMPACT

DESCRIPTION :



- Steel frames
- Laser-cut and bent sheet metal, tubes and sections
- Baked epoxy paint
- Ball-bearing wheels
- Total weight: 32 kg
- Maximum load: 150 kg
- Product life: 8 years* (excluding electrical part) *Subject to regular maintenance
- Lithium battery: 2 years warranty
- Lead battery: 1 year warranty

- Medical device (Regulation (EU) 2017/745)
- Electromagnetic compatibility
- Compliance with the following standards :
- NF EN 60601-1 : A1/2013
- NF EN 60601-1-2 : 2015
- NF EN ISO 10535 : 2021
- 24V low-voltage DC motor
- Power 24 V/120 VA
- Maximum thrust force: 3,000 N
- Electronic overload protection
- Protection class: II
- Type of protection: IP X4
- Contactor limit switches
- ABS housing
- Spiral cable remote control, low voltage 24V
- Overload protection by thermal cut-outs
- Sound level: 60 dB at 1 metre

MOTOR TECHNICAL FEATURES :

- Control box with built-in battery
- Battery life indicator light
- Power cord
- Emergency stop button
- 2-function remote control IPX4 protection
- Emergency lowering: electric and manual
- Mains cable connection for battery charging

DIMENSIONS :

- A.** Total length: 84 cm
- B.** Minimum useful height: 62 cm
- C.** Maximum useful height: 157 cm
- D.** Total height: 115 cm
- E.** Frame height: 8 cm
- F.** Minimum base width: 39 cm
- G.** Maximum base width: 60.5 cm
- H.** Turning diameter: 89 cm

