

BLUE EASYLEV 3 STAND UP LIFT

A budget-conscious device designed for high-end transfers!

The BLUE EASYLEV 3 **allows a semi-passive transfer.**

The patient must be able to put both feet flat.

It uses a patented movement that **takes care of the partial or complete lack of tonicity of the patient's legs and bust.**

The BLUE EASYLEV 3 supports **taller and/or lower toned patients** due to its increased range of motion and patented principle. The sling positioned at the pelvic level allows the patient to stand up by mechanically unfolding the lower body. The combined movement of the mast and the traction sling allows the patient to move from a sitting to a standing position while respecting their comfort and pathology.

It can also perform a subaxillary lift.

The quality of rolling has been improved thanks to its castors diameter 100 with double ball bearings.

The devices in the new range are equipped with two emergency lowering systems, one manual and one electric.

The positioning of the patient during the transfer allows access to the toilet, bathroom, etc.

This type of transfer allows the patient's psychological stimulation to continue, unlike a patient lift.

Lithium battery as standard to avoid additional operating costs.



WARRANTY
5
YEARS

CE



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE