

# BLUE WAYUP 3 STAND UP LIFT

## The ideal compromise between efficiency and cost!

The BLUE WAYUP 3 uses a patented movement that **takes care of the partial or complete lack of tonicity of the patient's legs.**

## The patient must be able to put both feet flat.

The sling is positioned at waist level, the adjustment of the hooks is done by putting the sling under tension with the patient's back against the backrest.

The quick and complete movement of the lifting arm **allows the patient to move from the sitting to the standing position, respecting the patient's comfort and pathology.**

Once the legs have been secured by the ergonomic knee support, the lifting movement can be performed. **The rolling quality has been improved** thanks to its castors diameter 100 with double ball bearings. The devices in the new range are equipped with two emergency lowering systems, one manual and one electric.

The positions thus obtained allow the activation of the patient's mobility and transit, facilitating care such as changing, washing and going to the toilet.

**Lithium battery as standard to avoid additional operating costs.**



NAUSICAA Medical provides a wide range of hoists slings for all transfers and compatible with all devices on the market.

WARRANTY  
**5**  
YEARS

CE



CONFORT | STYLE | BIEN-ÊTRE  
FABRICATION FRANÇAISE