

EASYLEV 5 STAND UP LIFT

A patented movement for unparalleled comfort!

The EASYLEV 5 uses a movement that takes care of the partial or complete lack of tonicity of the patient's legs and chest.

- **4 possible adjustments** : up, down, crossed (inner or outer handles).
- **3 different lengths + specific asymmetric positioning** (for hemiplegic patients for example).

The combined extension of the mast and the traction sling allows the patient to move from a sitting to a standing position while respecting their comfort and pathology.

The patient must be able to put both feet flat.

Once the legs are secured by the **ergonomic and adjustable knee support**, the lifting movement can be performed.

The design reduces the pressure exerted during the lifting phase.

The positions thus obtained allow the activation of the patient's mobility and transit, facilitating care such as changing, washing and going to the toilet.

Lithium battery as standard to avoid additional operating costs.



NAUSICAA Medical provides a wide range of hoists slings for all transfers and compatible with all devices on the market.

EASYLEV RANGE

SCAN TO WATCH



WARRANTY
5
YEARS



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE