

# NAUSIFLY 3

The NAUSIFLY range is designed for use in homes and institutions.

**It allows bed/chair transfers in a seated position** thanks to its 4-point parabolic spreader bar which increases patient comfort.

**It allows to pick up the patient from the ground.**

The devices in this range are equipped with two emergency lowering systems, one manual and one electric.

They have a wide mechanical or electrical opening legs and a folding frame.

**The rolling quality is optimal on all types of floor surfaces.**

**Lithium battery as standard to avoid additional operating costs.**

## Option



WEIGHT INDICATOR ACCESSORY



NAUSICAA Medical provides a wide range of hoists slings for all transfers and compatible with all devices on the market.



NEW

WARRANTY 5 YEARS



CONFORT STYLE BIEN-ÊTRE  
FABRICATION FRANÇAISE