

WAYUP 5 STAND UP LIFT

The most compact and quickest stand-up lift on the market!

The WAYUP 5 uses a patented movement that **takes care of the partial or complete lack of tonicity of the patient's legs.**

- **4 possible adjustments:** up, down, crossed (inner or outer handles).
- **3 different lengths + specific asymmetric positioning** (for hemiplegic patients for example).

The rapid and complete unfolding of the lifting arm allows the patient to move from the sitting to the standing position while respecting their comfort and pathology.

The patient must be able to put both feet flat.

Once the legs are secured thanks to its new adjustable and ergonomic knee support and its removable heel support, the lifting movement is possible.

The design reduces the pressure exerted during the lifting phase.

The resulting positions allow the activation of the patient's mobility and transit, facilitating care such as diapering, toileting and WCing.

Lithium battery as standard to avoid additional operating costs.



WARRANTY
5
YEARS



CONFORT STYLE BIEN-ÊTRE
FABRICATION FRANÇAISE