


NAUSIFLY 3

COMPACT

Specially designed **for use in small and confined spaces**. Easy to use.

Its advantages are:

- its **rolling quality**
- its **4-point comfort spreader bar**
- its **compact base**



Dimensions

- A. Total length: 104,5 cm
- B. Minimum useful height: 59 cm
- C. Maximum useful height: 164,5 cm
- D. Total height: 126 cm
- E. Frame height: 115 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 81 cm
- H. Turning diameter: 112 cm

REFERENCES

NAFLY3CPT-CBLI
NAFLY3CPT+-CBLI

NAUSIFLY 3


STANDARD

Its proportions make it particularly **easy to handle**.

Easy to use.

Its advantages are:

- its **rolling quality**
- its **4-point spreader bar**
- its **very large clearance**



Dimensions

- A. Total length: 124,5 cm
- B. Minimum useful height: 57,5 cm
- C. Maximum useful height: 179 cm
- D. Total height: 134,5 cm
- E. Frame height: 115 cm
- F. Minimum base width: 66 cm
- G. Maximum base width mechanical: 96 cm
Maximum base width electrical: 99 cm
- H. Turning diameter: 133,5 cm

REFERENCES

NAFLY3-CBLI
NAFLY3+-CBLI