



[www.nausicaa-medical.com](http://www.nausicaa-medical.com)

Before each use, it is important to check: the presence of the label, that the device does not show any tears, snags, holes or signs of fraying and the good condition of the hook and loop.

## HEEL PADS



### HEEL PADS SMALL MODEL


HELP TO PREVENT  
PRESSURE SORES  
ON HEELS

### HEEL PADS LARGE MODEL

HELP TO PREVENT  
PRESSURE SORES  
ON HEELS

#### Description :

The Heel Pads allow to to relieve pressure points and to wrap the foot in a furtrim, preventing rubbings on the heels and malleoli.



LMTAL-PF



The use of the device requires the presence and vigilance of a third person. When used by an agitated or disoriented person, monitoring becomes imperative.

This medical device must be implemented or trained by a health professional.



v 05.02/23

### Important Recommendations :

- In order to get a maximum efficiency of these sort of product, it is necessary :
- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient
- These products must not be in direct contact with a wounded skin.



## HELL PADS

COMBINATION OF POLYURETHAN FOAM  
AND JERSEY/COTTON POLYESTER FABRIC  
JERSEY 100% POLYESTER FABRIC  
100% POLYESTER FUR  
100% POLYESTER VELCRO STRAP

- Put the heel, fur side, right in the bottom of the heelpad.
- Close the heelpad on the foot through the velcro strap(s).