


WAYUP 5

COMPACT

Its very low ground clearance (6.5 cm) allows it to fit under **very low chairs, especially 6-wheel electric chairs.**
Its small footprint allows it to fit into any space (even confined ones, such as a toilet).



Dimensions

- A. Total length: 92 cm
- B. Minimum useful height: 74,5 cm
- C. Maximum useful height: 121 cm
- D. Total height: 100 cm
- E. Frame height: 6,5 cm
- F. Minimum base width: 34 cm
- G. Maximum base width: 56,5 cm
- H. Turning diameter: 90 cm


REFERENCES

WAYUP5-PF
 WAYUP5-LI-PF

WAYUP 5

MECHANICAL OR ELECTRICAL OPENING LEGS

Its wide base with mechanical or electrical spread allows an opening width of 90 cm making it **compatible with a large number of chairs.**
 Its removable platform makes it easy to use as a walking aid.



Dimensions

- A. Total length: 107 cm
- B. Minimum useful height: 82 cm
- C. Maximum useful height: 135 cm
- D. Total height: 106 cm
- E. Frame height: 11,5 cm
- F. Minimum base width: 54 cm
- G. Maximum base width mechanical: 83 cm
 Maximum base width electrical: 80 cm
- H. Turning diameter: 110 cm

REFERENCES

WAYUP5ECP-M-PF
 WAYUP5ECP-M-LI-PF
 W5ECP-E2-PF / W5ECP-E2-LI-PF