

## PRODUCT SHEET

### STAND ASSIST **EASYLEV**

### **COMPACT**

#### Description:



- Offset patient handlebar for optimal patient lifting comfort
- Adapts to all spaces and allows easy access to all types of areas thanks to its short base
- Fast and safe lifting in about thirty seconds
- Foldable, separable, and can be stored upright
- Steel chassis
- Ball bearing wheels
- Laser-cut and bent sheets, tubes, and profiles
- Oven-cured epoxy paint
- Total weight compact version: 28 kg
- Weight of the upper part: 16.40 kg (with one axle and two knobs)
- Weight of the Compact base: 11.60 kg (with the other axle and two knobs)
- Maximum load: 160 kg
- Product lifespan: 8 years\* (excluding electrical components)
- \* Subject to regular maintenance
- Lithium battery: 2-year warranty


- Medical devices (Regulation (EU) 2017/745)
- Compliance with the following standards:
  - NF EN 60601-1: 2006 A1/2013
  - NF EN 60601-1-2: 2015
  - NF EN ISO 10535: 2021
- Electromagnetic compatibility
- Low voltage 24V DC motors
- Power 24 V/120 VA
- Maximum push force: 3,000 N
- Electronic protection in case of overload
- Protection class: II
- Protection type: IP X4
- End-of-stroke stops by contactor
- ABS housing
- Spiral cable remote control, low voltage 24V
- Overload protection by thermal cut-outs
- Noise level: 60 dB at 1 meter

#### Motorization technical specifications:

- Control box with removable high-capacity battery
- Graduated battery charge indicator
- Autonomy displayed by indicator light
- Emergency stop button
- 2-function remote control, IPX4 protection (up/down)
- Emergency lowering: electric and manual
- Mains cord connection for battery charging

#### Dimensional characteristics:

- A.** Total length: 840 mm
- B.** Minimum usable height: 630 mm
- C.** Maximum usable height: 1580 mm
- D.** Total height: 1138 mm
- E.** Chassis height: 80 mm
- F.** Minimum base width: 390 mm
- G.** Maximum base width: 605 mm
- H.** Turning diameter: 890 mm
- I.** Ground clearance: 45 mm



**When stored folded upright** (see photo opposite):

- A. Width: 580 mm
- B. Depth: 570 mm
- C. Height: 1040 mm



**Packaging:** 1200x600x600 mm (stackable up to three products)

- 800x1200 mm pallet: up to 3 products
- 1000x1200 mm pallet: up to 6 products